


# Additional Overlay Options

## Category: Delta Pediatric Overlays

### Delta Foam

| Item #                  | Description                           | Qty/case |
|-------------------------|---------------------------------------|----------|
| DELTA-C11               | Overlay 32"x75"x3.5"                  | 1        |
| DELTA-120               | Overlay 34"x75"x3.5"                  | 1        |
| <b>Pediatric sizes:</b> |                                       |          |
| DELTA-200               | Isolette Matt Overlay 13.5"x25"x1.75" | 1        |
| DELTA-202               | Warmer Overlay 29"x23.75"x2.25"       | 1        |
| DELTA-204               | Crib Overlay 26.5"x53"x2.75"          | 1        |



Outstanding specialty surface for weights below 70 lbs. Ideal for pediatric settings.

### QUICK FACTS

#### Delta Pediatric Overlays

##### Primary Use:

- Skin management in infants and children

##### Key Advantages:

- Proven effectiveness\*
- High conformance to body contours
- Economical
- Protection in home or hospital


\*Studies available on request

## Category: Preventative Overlays

### Pressure Kair™

Channel cut side-to-side, then sliced head to foot to create individual support cells. Divided into distinctive head, torso, and foot zones.


| Item #    | Description                                | Qty/case |
|-----------|--|----------|
| SP789-000 | Pressure Kair 34"x72"x3.5" with Soft Skin® | 4        |
| SP889-000 | Pressure Kair 34"x80"x3.5" with Soft Skin® | 4        |



### PRT Matt™

Cross-sliced in head, torso, and foot zones. Distinctive lavender color.

| Item #    | Description  | Qty/case |
|-----------|--|----------|
| 50960-600 | PRT-MATT 34"x72"x3" cross cut with Soft Skin® (fitted) | 4        |
| 50960-604 | PRT-MATT 34"x72"x3" cross cut without Soft Skin®       | 4        |




### Span+Guard®

Comfortable, economical overlay. 3" peak, 2" base, in multi-zone design.

| Item #    | Description                                      | Qty/case |
|-----------|--|----------|
| SP799-000 | Span+Guard 6 zoned pad, 3"P, 2"B with Soft Skin® | 4        |

(Note: Covered under Medicare code E0199)



### QUICK FACTS

#### Preventative Overlays

##### Primary Use:

- Pressure ulcer prevention in low-mod risk patients


##### Key Advantages:

- Low cost
- Comfortable
- More conformance to body contours than innerspring mattresses

## Category: Convoluted Overlays

### Convoluted Overlays

| Item #    | Description                          | Qty/case |
|-----------|--------------------------------------|----------|
| SP22R-000 | 2" Convoluted Pads, Premium Density  | 12       |
| SP22S-000 | 2" Convoluted Pads, Standard Density | 12       |
| SP32R-000 | 3" Convoluted Pads, Premium Density  | 6        |
| SP35S-000 | 3" Convoluted Pads, Standard Density | 6        |
| SP45R-000 | 4" Convoluted Pads, Premium Density  | 4        |
| SP45S-000 | 4" Convoluted Pads, Standard Density | 4        |



Made from high-quality, combustion modified poly-urethane foam. Specifications vary to address any level of surface discomfort. All pads are 72"L x 34"W, and are compression rolled for efficient storage.

### QUICK FACTS

#### Convoluted Overlays

##### Primary Use:

- Comfort

##### Caution:

- Convoluted, or "eggcrate", overlays are recognized for comfort, not pressure management. Call your Span-America representative when selecting the pressure management support surface most effective for your patients' needs.

