


Encore Structural Highlights



Optional Accessories:

- Bed extenders to 84"
- Ergonomic 3-position assist device
- Bed & deck angle indicator
- Snap-in pendant holder fits all rails
- Molded plastic half-rail
- Headboard mount pendant holder
- Transport dollies
- Portable backup battery
- Regular backup battery
- Drainage hooks
- IV pole
- Trapeze

Height Floor to Mattress Support	7-7/8" - 31"
Deck Angles - Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 30°
Bed Length - Overall Length	82", 86", 90"
Bed Length - Retraction Travel	12"
Bed Length - Sleep Surface	76", 80", 84"
Bed Width - With Side Rails	40-1/2", 43-1/2", 46-1/2"
Bed Width - Sleep Surface	35", 39", 42"
Weight of Bed (without accessories)	253 lbs / 115 kg
Safe Working Load	600 lbs / 272 kg

Item #	Description
MC9R	Encore bed with GlideAlign retractability
MC9RW	Encore bed with GlideAlign retractability, with readywide expanders up to 42" wide
MC9RZ	Encore bed with GlideAlign retractability, with built-in battery backup
MC9RWZ	Encore bed with GlideAlign retractability, with readywide expanders up to 42" wide and built-in battery backup
R9955	Three position pivot assist - universal
R9956	Three position pivot assist - set
R9307H	Half head end side assist/rail - set
R9307H-L	Half head end side assist/rail - left
R9307H-R	Half head end side assist/rail - right
Q6691	Pivot assist wide mounting kit - 39"
Q6693	Pivot assist wide mounting kit - 42"
Q6690	Extender kit 4" - to 84"
Q6699	Dolly set
Q6698	Trapeze assembly
Q6705	Trapeze mounting frame
Q6708	Roller bumper kit
Q6703	Deck angle indicators with brackets

