

Seating

Category: Prevention without Positioning

Gel-T®

Segmented gel-filled bladder cradled between two soft foam layers. Provides excellent pressure redistribution and pain relief for those at moderate risk of pressure ulcers. Wipe-clean cover keeps out moisture and helps reduce shear and friction.

Cover also features non-skid bottom. Perfect for use in home or facility;

18-month warranty.

Height 2.5".

Weight limit 250 lbs.

* sizes 18" x 18" and larger

up to 350 lbs. (500 lbs. on Bariatric sizes)

Medicare Code:

General Use – E2601.

For sizes 22" or wider, use E2602.



Item #	Description	Qty/Cs
8015-05	Gel-T, 15"W x 15"L	1
8016-05	Gel-T, 16"W x 16"L	1
801618-05	Gel-T, 16"W x 18"L	1
8017-05	Gel-T, 17"W x 17"L	1
801816-05	Gel-T, 18"W x 16"L	1
8018-05	Gel-T, 18"W x 18"L	1
8019-05	Gel-T, 19"W x 19"L (fits Geri-chairs)	1
802018-05	Gel-T, 20"W x 18"L (fits Geri-chairs)	1
8021-05	Gel-T, 21"W x 21"L (fits Geri-chairs)	1
Bariatric sizes:		
B802218-05	Gel-T, 22"W x 18"L	1
B802418-05	Gel-T, 24"W x 18"L	1

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-05" suffix with "-29" in item number.

QUICK FACTS

Use for:

- Comfort
- Prevention of breakdown

RISK LEVEL

