

Geo-Mattress® Atlas®



Why "Progressive Resistance"?

Our design puts the most patient-friendly surface directly next to the skin, and the firmest, most dense material farthest away. The layers deflect in stages, allowing the Atlas to accommodate the widest possible range of weights and still remain comfortable.



The key features of the clinically proven Geo-Mattress series, enhanced to meet the unique pressure management challenges of the bariatric client. On top, the Geo-Matt design. Its 800 independent support cells protect against the damaging effects of friction and shearing – two major risk factors that increase dramatically when bed mobility and safe transfers are impeded by extreme weight. Beneath it, three high-resiliency layers offer progressively firmer levels of support. And at the foot, the proven Heel Slope® taper provides added protection for the heel and ankle. All protected by a fluid-proof cover and non-prorated 3-year warranty.

Additional bariatric options available – see page 18.

***Allow up to 10 working days lead time for these sizes.**

Weight limit: 600 lbs. (35"); 675 lbs. (39"); 750 lbs. (42", 48" and 54")

Warranty: 3 years, non-prorated **Mattress Weight:** 26 lbs. (35")

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

Description

✓ Bariatric Therapeutic Foam Mattress

QUICK FACTS

Indications:

- Prevention for high risk patients
- Patients 350-600 lbs

Key Advantages:

- Progressive resistance accommodates wide range of patient weights
- Geo-Matt® anti-shearing surface
- Heel Slope®

Item #	Description
A7535-29	Geo-Mattress™ Atlas, 75"L x 35"W x 7"H
A8035-29	Geo-Mattress™ Atlas, 80"L x 35"W x 7"H
A8435-29	Geo-Mattress™ Atlas, 84"L x 35"W x 7"H
* A8039-29	Geo-Mattress™ Atlas, 80"L x 39"W x 7"H
A8042-29	Geo-Mattress™ Atlas, 80"L x 42"W x 7"H
A8048-29	Geo-Mattress™ Atlas, 80"L x 48"W x 7"H
A8054-29	Geo-Mattress™ Atlas, 80"L x 54"W x 7"H
* A8442-29	Geo-Mattress™ Atlas, 84"L x 42"W x 7"H
* A8448-29	Geo-Mattress™ Atlas, 84"L x 48"W x 7"H
* A8454-29	Geo-Mattress™ Atlas, 84"L x 54"W x 7"H



SURFACES