

Risk Manager™

Category: Bedside Safety Mat

Minimize the Impact of Falls – On Your Residents and On Your Facility.

Statistics don't lie. Despite the best efforts of the healthcare community, falls by nursing home residents and hospital patients do occur. When they do, the physical, legal, and financial impact can be devastating to all involved.


Span-America's new Risk Manager™ bedside safety mat is designed to help.

The Risk Manager is engineered from an ultra-tough, impact-absorbing lightweight elastomeric gel compound commonly used in expensive orthotics and premium footwear.


Independent laboratory testing (see chart) verifies that the Risk Manager greatly reduces the peak G-forces generated by falls, easily surpassing the impact reduction provided by thicker, bulkier gym mat-style pads. It also outperforms the most commonly purchased specialty fall strip, which relies on the same low-tech foam scrap padding typically found under residential carpeting.

IMPACT TEST COMPARISON

Drop Height 36"



¹Source: National Technical Systems, Detroit, Michigan. 2008. Results on file.



Premium elastomeric gel surface reduces impact from falls. Bonus: Provides welcome comfort for caregivers.



Nano-silver throughout protects mat from stains, odor-causing bacteria, mold and mildew. Mat wipes clean with most standard cleaning agents.



Slide-resistant bottom surface grips floor, helps keep mat in place. **Textured top** provides added safety.



Low profile and beveled edge allow easy on and off for wheelchairs and lifts.



Not a trip hazard, unlike bulky gym-style mats. Thin, stable surface means staff won't move Risk Manager from beside bed.

Item #	Risk Manager™ Bedside Safety Mat	Qty/Cs
37052	24"W x 70"L x 1"H, standard – single	1/case
94612	24"W x 70"L x 1"H, standard – bulk pack	5/case
8705236	36"W x 70"L x 1"H, XL – single	1/case
9611836	36"W x 70"L x 1"H, XL – bulk pack	5/case

Warranty: 1 year
Weight: 6 lbs. (24"), 9.5 lbs. (36")
Flammability: Fire retardant – passes Cal TB 117

Stain- and oil-resistant
Fluidproof, will not absorb moisture
Non-toxic. Latex-free. Lead-free.
Rolls up for easy transport, storage and shipping.

Optional Accessory*

RISK MANAGER™ GRIP STRIPS

For use on highly-polished floors where extra slide resistance may be desired.

*(Included as standard with XL version.)

Item # 43398 (3/Package)
 (2" wide x 14" long)



Now Available!
XL Version
36" Wide



SURFACES



SEATING



POSITIONERS



SKIN CARE

Risk Manager

Bedside Safety Mat

5-Way Protection:

- **Impact-reducing**, proprietary elastomeric gel
- **Textured top** helps reduce slipping
- **Slide-resistant** bottom surface
- **Fire-resistant**, meets Cal TB 117

PLUS:

Infused with Nano-silver

Independently documented, superior impact reduction provides an extra measure of bedside protection for fragile, at-risk residents.

SPANAmerica
Innovative Solutions.

Minimize the impact of falls –

On your residents *and* on your facility.

Statistics don't lie. Despite the best efforts of the healthcare community, falls by nursing home residents and hospital patients do occur. When they do, the physical, legal, and financial impact can be devastating to all involved.

Span-America's new Risk Manager™ bedside safety mat is designed to help.

Science. Safety. Simplicity.


The Risk Manager is engineered from an ultra-tough, impact-absorbing lightweight elastomeric gel compound commonly used in expensive orthotics and premium footwear.

Independent laboratory testing¹ (see chart) verifies that the Risk Manager greatly reduces the peak G-forces generated by falls, easily surpassing the impact reduction provided by thicker, bulkier gym mat-style pads. It also outperforms the most commonly purchased specialty fall strip, which relies on the same low-tech foam scrap padding typically found under residential carpeting.

Best of all, the Risk Manager provides this performance while also delivering crucial advantages in cleanability, caregiver comfort, ease of use, and portability.

IMPACT TEST COMPARISON

Drop Height 36"



¹Source: National Technical Systems, Detroit, Michigan. 2008. Results on file.




Premium elastomeric gel surface reduces impact from falls. Bonus: Provides welcome comfort for bedside caregivers.



Nano-silver throughout. Protects mat from stains, odor-causing bacteria, mold and mildew. Mat wipes clean with most standard cleaning agents.



Slide-resistant bottom grips floor, helps keep mat in place. **Textured top** provides added safety.



Low profile and beveled edge allow easy on and off for wheelchairs and lifts.



Not a trip hazard, unlike bulky gym-style mats. Thin, stable surface means staff won't move Risk Manager from beside bed.

Optional Accessory* RISK MANAGER™ GRIP STRIPS

For use on highly-polished floors where extra slide resistance may be desired



*(Included as standard with XL version.)

Item # 43398 (3/Case)
(2" wide x 14" long)

Item #	Risk Manager Bedside Safety Mat	Qty/Cs
37052	24"W x 70"L x 1"H, standard – single	1/case
94612	24"W x 70"L x 1"H, standard – bulk pack	5/case
8705236	36"W x 70"L x 1"H, XL – single	1/case
9611836	36"W x 70"L x 1"H, XL – bulk pack	5/case
Warranty:	1 year	
Weight:	6 lbs. (24"), 9.5 lbs. (36")	
Flammability:	Fire retardant – passes Cal TB 117	
	Stain- and oil-resistant	
	Fluid-resistant	
	Non-toxic. Latex-free. Lead-free.	
	Rolls up for easy transport, storage and shipping.	