

Soft-Touch® Sitters

Supportive, Stylish Seating

Base Styles

Three styles to fit your needs



Sitter with Stationary Base



Sitter with Mobile Base



Floor Sitter Kit for Sizes 1,2,3

Features



Sleek modern design provides up to 25° of tilt




Attachment straps for securing to a standard chair



Adjustable 8-point harness



Contoured head and lateral support, shown as top view



Built-in abductor



Optional Headrest

Soft-Touch[®] Sitters

5 year warranty

Sitter Only

5 sizes available for that “custom” fit



Size 1



Size 2



Size 3



Size 4



Size 5

Applications

Sitters can be used in a variety of ways



Attach to a chair



Jogger for Sizes 1,2 (page 22)



EIO Push Chair for Sizes 1,2 (page 20)

Sitter Specifications	Size 1	Size 2	Size 3	Size 4	Size 5
Target Weight	20 - 30 lbs.	30 - 60 lbs.	40 - 80 lbs.	50 - 110 lbs.	90 - 200 lbs.
Target Height	30" - 35"	35" - 48"	40" - 56"	50" - 65"	54" - 74"
Target Age	1 - 2 years	2 - 6 years	6 - 9 years	9 - 14 years	13 years - Adult
Inside Seat Depth	7"	9.5"	11.5"	13.5"	16"
Inside Seat Width	8"	10"	12"	14"	16"
Seat Back Height	21"	26"	29"	35"	39"
Floor to Seat Ht. UPRIGHT on Stationary Base	10"	10"	9.75"	11.5"	12.5"
Floor to Seat Ht. TILTED on Stationary Base	10.5"	10.5"	10.5"	13"	14.25"
Floor to Seat Ht. UPRIGHT on Mobile Base	12.5"	13.25"	13"	14.75"	16"
Floor to Seat Ht. TILTED on Mobile Base	12.75"	13.75"	13.75"	16"	17.75"
Floor to Seat Ht. UPRIGHT on Floor Base	5"	5.5"	5.5"	N/A	N/A
Floor to Seat Ht. TILTED on Floor Base	4.5"	4.5"	5"	N/A	N/A
Overall Seat Depth	12"	14"	16"	19"	21"
Overall Seat Width	12"	14"	16"	18"	20"
Overall Seat Height	25"	29"	34"	40"	44"
Weight of Seat	6 lbs.	10 lbs.	13 lbs.	15 lbs.	19 lbs.

Base Specifications	Floor Base	Stationary (Small)	Stationary (Large)	Mobile Base (Small)	Mobile Base (Large)
Use with Sitter Sizes	1, 2, 3	1, 2, 3	4, 5	1, 2, 3	4, 5
Overall Width	19"	17.5"	19.25"	17.5"	19.25"
Overall Depth	22.5"	30"	35"	30"	35"
Overall Height	12.5"	16"	23"	19.5"	26.5"
Weight of Base	14 lbs.	12 lbs.	15 lbs.	14 lbs.	17 lbs.
Degree of Tilt	25°	25°	25°	25°	25°

Soft-Touch® Sitters (w/ Gray Harness)	Lilac	Aqua	Dark Gray
Size 1 Sitter Only	77010102	77010104	77010109
Size 2 Sitter Only	77010202	77010204	77010209
Size 3 Sitter Only	77010302	77010304	77010309
Size 4 Sitter Only	77010402	77010404	77010409
Size 5 Sitter Only	77010502	77010504	77010509
Size 1 Sitter on Stationary Base	77110102	77110104	77110109
Size 2 Sitter on Stationary Base	77110202	77110204	77110209
Size 3 Sitter on Stationary Base	77110302	77110304	77110309
Size 4 Sitter on Stationary Base	77110402	77110404	77110409
Size 5 Sitter on Stationary Base	77110502	77110504	77110509
Size 1 Sitter on Mobile Base	77210102	77210104	77210109
Size 2 Sitter on Mobile Base	77210202	77210204	77210209
Size 3 Sitter on Mobile Base	77210302	77210304	77210309
Size 4 Sitter on Mobile Base	77210402	77210404	77210409
Size 5 Sitter on Mobile Base	77210502	77210504	77210509
Size 1 Floor Sitter	78010102	78010104	78010105
Size 2 Floor Sitter	78010202	78010204	78010209
Size 3 Floor Sitter	78010302	78010304	78010309

Optional Items	Part No.
Optional Floor Sitter Wedge (fits sizes 1-3 sitter)	78000009
Optional Small Stationary Base (fits sizes 1-3 sitter)	77100100
Optional Small Mobile Base (fits sizes 1-3 sitter)	77200100
Optional Headrest	77020000
Optional Large Stationary Base (fits sizes 4-5 sitter)	77100200
Optional Large Mobile Base (fits sizes 4-5 sitter)	77200200