


# MPS - Customize Your Size

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!  
Please mark clearly with an X the items you wish to purchase.

Account #: \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

## Special Tomato® MPS Features:

- Individualized adjustment allows for many years of growth
- Seamless latex-free cushions are impermeable to fluids and soft, yet durable
- Cleans easily with common disinfectants
- All Special Tomato MPS Seating Systems are FMVSS213 (Crash Testing) Compliant



## Select Your Cushions: Choose One of Each

MPS Seat Specifications	Small Seat Cushions				Large Seat Cushions			
	MS1	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	9 mos - 3 yrs	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	28" - 37"	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Weight	80 lbs.				130 lbs.			
Seat Width	7"	9"	10"	10"	11"	12"	13"	13"
Seat Depth	6"	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16"				8" - 18"			
Hip Flex Angle	-15° to +15°				-15° to +15°			
Floor to Seat Height	15"				18"			

MPS Back Specifications	Small Back Cushions		Large Back Cushions		
	MS2	MS4	ML2	ML4	ML6
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult	10 yrs - Adult
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult	48" - Adult
Maximum User Weight	80 lbs.		150 lbs.		
Back Width	8"	9"	10"	11"	13"
Back Height	8.5"	11"	11"	13"	13"
Top of Back Cushion	11" - 16"	15" - 20"	15" - 20"	19" - 22"	19" - 22"
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"	26" - 30"

Product	Small Dimensions	Large Dimensions	Small Weight	Large Weight
Tray	24"W × 17"D × 2"H		8 lbs.	
MPS Seat Only	15"W × 19"D × 42"H	19"W × 23"D × 47"H	21 lbs.	26 lbs.
MPS Floor Kit	19"W × 37"D × 37"H	22"W × 41"D × 42"H	39 lbs.	44 lbs.
MPS Push Chair Kit	25"W × 38"D × 42"H	25"W × 42"D × 47"H	68 lbs.	77 lbs.

\* ML2 designed to fit both the Small and Large system, providing added wider Back Cushion option for Small size systems.

# MPS - Customize Your Size

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!  
Please mark clearly with an X the items you wish to purchase.

## Customize Your Size



Small Shell (MS)



Large Shell (ML)


### Shells

Choose a Shell size to determine your options

### Options

#### Headrests Cushions for both Small (MS) and Large Shells (ML)

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)



7"W x 6"H

MBH (Basic)



7.5"W x 6.5"H  
1.5"D

MBL (Basic w/ Laterals)



7"W x 6"H  
2.5"D

#### Back Cushions for Small Shell (MS)



8"W x 8.5"H




9"W x 11"H



10"W x 11"H

#### Back Cushions for Large Shell (ML)



10"W x 11"H



11"W x 13"H




13"W x 13"H


#### Seat Cushions for Small Shell (MS)



7"W x 6"D



9"W x 8"D




10"W x 10"D




10"W x 12"D


#### Seat Cushions for Large Shell (ML)



11"W x 11"D



12"W x 13"D



13"W x 15"D



13"W x 17"D

# MPS - Customize Your Size

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!  
Please mark clearly with an X the items you wish to purchase.

## Customize Your Size



### MPS Seat Shell

Small - 79040109  
\$1,306.45  
 Large - 79040309  
\$2,388.49



### MPS Shell with Footrest

Small - 79030109  
\$2,893.12  
 Large - 79030309  
\$3,928.52



### MPS Mobile Floor Base

Small - 79010109  
\$4,873.28  
 Large - 79010309  
\$6,295.40

## Select MPS Cushions: One of Each

### Headrest Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MHS	<input type="checkbox"/> 79100051	<input type="checkbox"/> 79100052	<input type="checkbox"/> 79100054	<input type="checkbox"/> 79100059	\$236.29
MHL	<input type="checkbox"/> 79110051	<input type="checkbox"/> 79110052	<input type="checkbox"/> 79110054	<input type="checkbox"/> 79110059	\$294.42
MBH	<input type="checkbox"/> 79130051	<input type="checkbox"/> 79130052	<input type="checkbox"/> 79130054	<input type="checkbox"/> 79130059	\$236.29
MBL	<input type="checkbox"/> 79140051	<input type="checkbox"/> 79140052	<input type="checkbox"/> 79140054	<input type="checkbox"/> 79140059	\$294.42

### Back Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MS2	<input type="checkbox"/> 79120551	<input type="checkbox"/> 79120552	<input type="checkbox"/> 79120554	<input type="checkbox"/> 79120559	\$331.26
MS4	<input type="checkbox"/> 79120151	<input type="checkbox"/> 79120152	<input type="checkbox"/> 79120154	<input type="checkbox"/> 79120159	\$360.50
ML2	<input type="checkbox"/> 79120751	<input type="checkbox"/> 79120752	<input type="checkbox"/> 79120754	<input type="checkbox"/> 79120759	\$412.73
ML4	<input type="checkbox"/> 79120351	<input type="checkbox"/> 79120352	<input type="checkbox"/> 79120354	<input type="checkbox"/> 79120359	\$423.81
ML6	<input type="checkbox"/> 79120451	<input type="checkbox"/> 79120452	<input type="checkbox"/> 79120454	<input type="checkbox"/> 79120459	\$449.24

\* ML2 designed to fit both the Small and Large system, providing added wider Back Cushion option for Small size systems.

### Seat Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MS1	<input type="checkbox"/> 79050551	<input type="checkbox"/> 79050552	<input type="checkbox"/> 79050554	<input type="checkbox"/> 79050559	\$331.26
MS3	<input type="checkbox"/> 79050151	<input type="checkbox"/> 79050152	<input type="checkbox"/> 79050154	<input type="checkbox"/> 79050159	\$360.50
MS5	<input type="checkbox"/> 79050251	<input type="checkbox"/> 79050252	<input type="checkbox"/> 79050254	<input type="checkbox"/> 79050259	\$389.73
MS7	<input type="checkbox"/> 79050651	<input type="checkbox"/> 79050652	<input type="checkbox"/> 79050654	<input type="checkbox"/> 79050659	\$412.73
ML1	<input type="checkbox"/> 79050751	<input type="checkbox"/> 79050752	<input type="checkbox"/> 79050754	<input type="checkbox"/> 79050759	\$412.73
ML3	<input type="checkbox"/> 79050351	<input type="checkbox"/> 79050352	<input type="checkbox"/> 79050354	<input type="checkbox"/> 79050359	\$423.81
ML5	<input type="checkbox"/> 79050451	<input type="checkbox"/> 79050452	<input type="checkbox"/> 79050454	<input type="checkbox"/> 79050459	\$434.88
ML7	<input type="checkbox"/> 79050851	<input type="checkbox"/> 79050852	<input type="checkbox"/> 79050854	<input type="checkbox"/> 79050859	\$460.55

\* MS1 Cushion is not certified for use in a motor vehicle.

\$

Your MSRP Total

# MPS - Customize Your Size

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!  
Please mark clearly with an X the items you wish to purchase.

## MPS Accessories



### Mobile Floor Base

Built-in tray storage on back. Four Swivel casters with rear locks. 10° - 25° of tilt.

Size	Part #	MSRP
Small	___ 79070100	\$2,063.75
Large	___ 79070300	\$2,641.21



### MPS Tray

Depth adjustable. 60° tilt adjustment. Tray easily attaches to back of Mobile Floor Base for storage.

Size	Part #	MSRP
Small	___ 79060159	\$643.00
Large	___ 79060359	\$725.13