

Special Tomato CatchUp Chair™ Order Form



The seat that “catches up” with your child’s growth.

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!

Please mark clearly with an X the items you wish to purchase.



Account #: _____ Ship to: _____
 Date: _____ Name: _____
 PO #: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

CatchUp Chair Specifications

ITEM	ENGLISH	METRIC
User Weight Capacity	250 lbs	91 kg
Armrest Pad- SMALL	10"	25.4 cm
Seat Width (Between Armrests) SMALL	14"	35.5 cm
Armrest Pad- LARGE	14"	35.5 cm
Seat Width (Between Armrests) LARGE	15.5"	39.3 cm
Seat Width (No Armrests)	16"	40.6 cm
Seat Depth	11" to 17"	27.9 cm to 43.2 cm
Seat Back Height	*13.5" to variable	34.3 cm
Backrest Extension Height (off the seat cushion)	14" to 28"	35.5 cm to 71.1 cm
Winged Headrest Height Range	10" to 28"	25.4 cm to 71.1 cm
Shoulder Strap Height	13.5" to variable	33 cm
Seat to Footrest (variable)	4.5" to 18"	11.4 cm to 45.7 cm
Seat to Floor	13" to 29"	33 cm to 73.7 cm
Seat to Back Angle	100°	
Seat Angle to Floor	5°	
Chair Weight	75 lbs	34 kg
Push Handle Height	34"	86.4 cm
Overall Width	29"	73.7 cm
Overall Height (Top of Seatback to Floor)	42"	106.7 cm
Overall Length	37"	94 cm
Caster Size	5"	12.7 cm

#60000000 shown here with Backrest Extension (sold separately)

Adjustable 8-point chest harness (large & small Included)

Shown here with Tray & Arm Rests (sold separately)



No-Tool Seat & Foot Plate Height & Depth Adjustment

Four locking casters




PART #	DESCRIPTION	MSRP	QTY
60000000	CatchUp Chair Includes: 1" Padded Seat and Back Cushion + (Lg & Sm 8 Point Harness)	\$3,748.45	
60030109	CatchUp Adjustable Armrests (Pair) SMALL (Required for Tray)	\$548.26	
60030209	CatchUp Adjustable Armrests (Pair) LARGE (Required for Tray)	\$646.09	
60020109	CatchUp Chair Tray Kit for SMALL Armrests (Small Armrests Required / sold separately)	\$548.26	
60020209	CatchUp Chair Tray Kit LARGE Armrests (Large Armrests Required / Sold Separately)	\$646.09	
60040000	CatchUp Backrest Extension	\$318.16	
60050009	CatchUp Winged Headrest (Requires Backrest Extension)	\$283.43	



60020109 / 60020209



60030109 / 60030209



60040000



60050009

Can be used with up to a size 4 Soft-Touch® Liners

Seat Cushions	Cherry	Lilac	Aqua	Dark Gray	Inner Dimensions	Outer Dimensions	MSRP	CHECK BOX TO ORDER
Size 1 (LB1) *	72050101	72050102	72050104	72050109	7.5"W x 10"H	10.5"W x 10"H	\$385.95	
Size 2 (LS2)	72040201	72040202	72040204	72040209	10"W x 10"D x 4"H	12"W x 13"D	\$419.50	
Size 3 (LS3)	72040301	72040302	72040304	72040309	12"W x 12"D x 4.5"H	14"W x 15"D	\$479.09	
Size 4 (LS4)	72040401	72040402	72040404	72040409	14" W x 14" D X 4.5" H	17.5" W x 17.5" D	\$550.96	

Back Cushions	Cherry	Lilac	Aqua	Dark Gray	Inner Dimensions	Outer Dimensions	MSRP	CHECK BOX TO ORDER
Size 1 (LB1) *	72050101	72050102	72050104	72050109	7.5"W x 10"H	10.5"W x 10"H	\$385.95	
Size 2 (LB2) *	72050201	72050202	72050204	72050209	8.5"W x 13.5"H	12.25"W x 13.5"H	\$419.50	
Size 3 (LB3) *	72050301	72050302	72050304	72050309	10"W x 14.5"H	13.25"W x 14.5"H	\$479.09	
Size 4 (LB4) *	72050401	72050402	72050404	72050409	13" W x 15.5" H	14.5"W x 15.5" H	\$550.96	

Can be used with up to a size 3 Soft-Touch® Sitter

Soft-Touch® Sitters (w/ Gray Harness)	Lilac	Aqua	Dark Gray	MSRP	CHECK BOX TO ORDER
Size 1 Sitter Only	77010102	77010104	77010109	\$814.74	
Size 2 Sitter Only	77010202	77010204	77010209	\$1,072.04	
Size 3 Sitter Only	77010302	77010304	77010309	\$1,363.67	

Sitter Specifications	Size 1	Size 2	Size 3
Target Weight	20 - 30 lbs.	30 - 60 lbs.	40 - 80 lbs.
Target Height	30" - 35"	35" - 48"	40" - 56"
Target Age	1 - 2 years	2 - 6 years	6 - 9 years
Inside Seat Depth	7"	9.5"	11.5"
Inside Seat Width	8"	10"	12"
Seat Back Height	21"	26"	29"
Overall Seat Depth	12"	14"	16"
Overall Seat Width	12"	14"	16"
Overall Seat Height	25"	29"	34"
Weight of Seat Overall Depth	6 lbs.	10 lbs.	13 lbs.



Shown with size 3 Soft-Touch® Sitter