

StanzMe PRONE Stander October 2022

Address: **16 Signal Hill Road, Brookfield, CT 06804** | Contact number: **1 (888) 714-8405** | Email: **sales@StanzMe.com**

Mark For: Date: Dealer Acct #: Dealer: Dealer Contact: Dealer Address: Dealer City: Dealer Phone: Email: Confirm Via:	Submitting for: Quote Order PO#: ADDITIONAL SHIPPING INFORMATION Ship To: Attention: Address: Address: Ship To City: ST: ZIP: Ship To Phone: Fax:
Fax Email	

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

StanzMe Prone Stander - Package

Package includes

1. Stander Frame

2. Tray

3. Padded Chest Board

4. Lateral Trunk Supports (Pair)

5. Lateral Hip Supports (Pair)

6. Adjustable Pelvic Support Pad

7. Lateral Knee Supports (Pair)

8. Shoe Holder

9. Two 2" Velcro Straps

10. Backrest Extension Pad (**optional - additional cost**)

Item Code	Description	HCPCS	MSRP
SPS-4000S	StanzMe PRONE Stander Package - Small	E0638	\$2,800.00
SPS-4000L	StanzMe Prone Stander Package - Large	E0638	\$3,000.00

StanzMe PRONE Stander - Itemized Pricing

StanzMe PRONE Stander - Small				
	Item Code	Description	HCPCS	MSRP
X	SPS-FRAME-SMALL	Frame U-Shaped Ergo Design	E0638	\$1000.00
X	SPS-UCESS	Tray + adjustable bracket + hardware	E1399	\$400.00
X	SPS-PCB	Padded Chest Board	E1399	\$200.00
X	SPS-LTS	Lateral Trunk Supports (Pair) + Hardware	E1399	\$200.00
X	SPS-LHS	Lateral Hip Supports (Pair) + Hardware	E1399	\$200.00
X	SPS-APSP	Adjustable Pelvic Support Pad + Hardware	E1399	\$400.00
X	SPS-LKS	Lateral Knee Supports (Pair) + Hardware	E1399	\$200.00
X	SPS-SH	Shoe Holder adjustable	E1399	\$100.00
X	SPS-VS	Ankle and Toe Straps 1" 2 sets	E1399	\$100.00
	SPS-BKSE	Back Support ext pad (extra charge) \$200 (not included in the pkg)	E1399	Optional charge
Total (shipping additional)				\$2,800.00
ALL PARTS ARE INCLUSIVE IN THE PACKAGE EXCEPT SPS-BKSE				

StanzMe PRONE Stander - Itemized Pricing

Stanzme Prone Stander - Large

	Item Code	Description	HCPCS	MSRP
X	SPS-FRAME-LARGE	Frame - U-shaped Ergo design	E0638	\$1200.00
X	SPS-UCESS	Tray + Adjustable angle + hardware	E1399	\$400.00
X	SPS-PCB	Padded Chest Board	E1399	\$200.00
X	SPS-LTS	Lateral Trunk Supports (Pair) + Hardware	E1399	\$200.00
X	SPS-LHS	Lateral Hip Supports (Pair) + Hardware	E1399	\$200.00
X	SPS-APSP	Adjustable Pelvic Support Pad + Hardware	E1399	\$400.00
X	SPS-LKS	Lateral Knee Supports (Pair) + Hardware	E1399	\$200.00
X	SPS-SH	Shoe Holder adjustable	E1399	\$100.00
X	SPS-VS	Ankle and Toe Straps 1" 2 sets	E1399	\$100.00
	SPS-BKSE	Back Support ext pad (extra charge) \$200 (not included in the pkg)	E1399	Optional charge
Total (shipping additional)				\$3,000.00

ALL PARTS ARE INCLUSIVE IN THE PACKAGE EXCEPT SPS-BKSE

Technical Specifications - Metric

Stanzme Prone Vertical Stander

Size	Small	Large
User age	9 months - 5 years	6 - 13 years
User height (approximate)	24" - 36"	34" - 48"
Footrest dimensions	22"W x 12"D/56cm x 32cm	25"W x 16"D/64cm x 40cm
Base - width x length (fully extended)	22" x 32"/56cm x 82cm	26" x 45"/66cm x 115cm
Footrest to top of Axilla	17" - 29"/56cm - 74cm	26" - 36"/66cm - 92cm
Tray size	17" x 24"/44cm x 61cm	19" x 28"/48cm x 70cm
Floor to top of chest board	27"/71cm	42"/107cm

Standing Measurements:

A ___ Knee height (back of knee to floor)

B ___ Hip height (top of hip bone to floor)

C ___ Elbow height (bent elbow to floor)

D ___ ***Axilla height** (armpit to floor)

E ___ Shoulder height (top of shoulder to floor)

F ___ Overall height (top of head to floor)


G ___ Hip width (widest location)

H ___ Trunk width (widest location)

***Key Measurement** –Should correspond to specification for footrest to top of chest board

Additional Accessories:

Backrest Extension Pad



Prone Stander Sizes



9 mo.- 5 years



6-12 year