

StanzMe Supine Stander October 2022

Address: **16 Signal Hill Road, Brookfield, CT 06804** | Contact number: **1 (888) 714-8405** | Email: **sales@StanzMe.com**

Mark For: Date: Dealer Acct #: Dealer: Dealer Contact: Dealer Address: Dealer City: Dealer Phone: Email: Confirm Via:	Submitting for: Quote Order PO#: ADDITIONAL SHIPPING INFORMATION Ship To: Attention: Address: Address: Ship To City: ST: ZIP: Ship To Phone: Fax:
Fax Email	

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

StanzMe Supine Stander - Package

Package includes

1. Stander Frame

2. Tray

3. Chest Harness

4. Lateral Trunk Supports (Pair)

5. Lateral Hip Supports (Pair)

6. Flat Headrest with Lateral Support

7. Lateral Knee Supports (Pair)

8. Shoe Holder

9. Two 2" Velcro Straps

10. Padded Back Board

Item Code	Description	HCPCS	MSRP
SSS-6000S	StanzMe Supine Stander Package - Small	E0638	\$3,200.00
SSS-6000L	StanzMe Supine Stander Package - Large	E0638	\$3,400.00

StanzMe SUPINE Stander - Itemized Pricing

StanzMe SUPINE Stander - Small				
	Item Code	Description	HCPCS	MSRP
X	SSS-FRAME-SMALL	Frame	E0638	\$1200.00
X	SSS-UESS	Tray with adjustable angle hardware	E1399	\$400.00
X	SSS-PBB	Padded Back Board	E1399	\$200.00
X	SSS-LTS	Lateral Trunk Supports (Pair)	E1399	\$200.00
X	SSS-LHS	Lateral Hip Supports (Pair)	E1399	\$200.00
X	SSS-HLHP	Headrest with Adjustable Lateral Supports	E1399	\$400.00
X	SSS-CH	Chest Harness	E1399	\$200.00
X	SSS-LKS	Lateral Knee Supports (Pair)	E1399	\$200.00
X	SSS-SH	Shoe Holder	E1399	\$100.00
X	SSS-VS	Two PAIR of 1" Velcro Straps	E1399	\$100.00
Total (shipping additional)				\$3,200.00

StanzMe Supine Stander - Large				
	Item Code	Description	HCPCS	MSRP
X	SSS-FRAME-LARGE	Frame	E0638	\$1,400.00
X	SSS-UESS	Tray with adjustable angle hardware	E1399	\$400.00
X	SSS-PBB	Padded Back Board	E1399	\$200.00
X	SSS-LTS	Lateral Trunk Supports (Pair)	E1399	\$200.00
X	SSS-LHS	Lateral Hip Supports (Pair)	E1399	\$200.00
X	SSS-HLHP	Headrest with Adjustable Lateral Supports	E1399	\$400.00
X	SSS-CH	Chest Harness Butterfly	E1399	\$200.00
X	SSS-LKS	Lateral Knee Supports (Pair)	E1399	\$200.00
X	SSS-SH	Shoe Holder	E1399	\$100.00
X	SSS-VS	Two PAIR of 2" Velcro Straps	E1399	\$100.00

Total **(shipping additional)** \$3,400.00

**ALL PARTS ARE INCLUSIVE IN THE PACKAGE EXCEPT Power actuator kit in LARGE unit ONLY
Power actuator conversion kit available for LARGE SUPINE at additional \$200 cost plus shipping**

Technical Specifications

StanzMe SUPINE Vertical Stander

Size	Small	Large
User age	1 - 5 years	7 - 12 years
User height (approximate)	33" - 55"	53" - 65"
Footrest dimensions	22"W x 12"D/56cm x 32cm	25"W x 16"D/64cm x 40cm
Base - width x length (fully extended)	22" x 42"/55cm x 105cm	29.5" x 45.5"/75cm x 116cm
Footrest to top of chest board	33.5" - 55"/85cm - 140cm	53" - 59"135cm - 150cm
Tray size	17" x 24"/44cm x 61cm	19" x 28"/48cm x 70cm
Range to Supine in Degrees	180 degrees to 90 degrees	180 degrees to 90 degrees


Standing Measurements:

- A ___ Knee height (back of knee to floor)
- B ___ Hip height (top of hip bone to floor)
- C ___ Elbow height (bent elbow to floor)
- D ___ * **Axilla height** (armpit to floor)
- E ___ Shoulder height (top of shoulder to floor)
- F ___ Overall height (top of head to floor)
- G ___ Hip width (widest location)
- H ___ Trunk width (widest location)

***Key Measurement** –Should correspond to specification for footrest to top of chest board

Additional Accessories:

- Power actuator conversion kit



Supine Stander Sizes



- 1 - 7 years (small)



- 7 - 14 years (large)