

# SIMPLY OUTSTANDING

Stanzme is designed for daily standing therapeutic needs of children from early intervention to age 14. Its simple structure and application is easy to use by parents who may not have enough experience with this type of equipment. It is simple in design and construction as compared to others and allows caregivers to position the child quickly and with ease.



# SUPINE STANDER

**ADDRESS**

16 Signal Hill Rd  
Brookfield, CT 06804

**TELEPHONE**

+(888) 714-8405

**EMAIL**

[sales@stanzme.com](mailto:sales@stanzme.com)

**WEB**

[www.stanzme.com](http://www.stanzme.com)



# Supine Stander

Stanzme Supine Stander is a posterior to vertical stander. It enables children with low muscle tone and poor head control to bear weight and stay upright against gravity. This stander is designed for easy supine transfer while assisting children to an upright position.

## Features

- Lateral trunk and hip pads
- Padded back board
- Trunk and hip straps
- Lateral knee supports with straps
- Shoe holder
- Tray with adjustable angle and removable bracket
- Locking casters
- Anterior chest harness
- Headrest with lateral adjustable pads





Item	SSS-6000S	SSS-6000L
User Age/Weight	1-7 years/ 70lbs	7-14 years/ 130lbs
User Height (approximate)	33" - 55" / 84cm - 140cm	53" - 65" / 135cm - 165cm
Footrest dimensions	22"W x 12"D / 56cm x 32cm	25"W x 16"D / 64cm x 40cm
Base-width x length (fully extended)	21.5" x 41.5" / 55cm x 105cm	29.5" x 45.5" / 75cm x 116cm
Footrest to top of supine back	33.5" x 55" / 85cm x 140cm	53" x 59" / 135cm x 150cm
Tray Size	17" x 24" / 44cm x 61cm	19" x 28" / 48cm x 70cm
Range to Supine	Supine to 90°	Supine to 90°