



# SHORT CONTOUR CUSHION


The Short Contour Cushion (2" cell height) is recommended for the active individual that can perform independent weight shifts, or for those individuals who are at low risk for pressure ulcer. The Short Contour is the perfect cushion for sports enthusiasts as well, providing the ideal seating surface for your everyday chair, sports chair or other sports equipment.



2" cell height  
HCPCS CODE - K0734



HCPCS CODE - K0736 Dual Valve



HCPCS CODE - K0736 Four Valve



2" cell height with cover

- 2" high cell cushion is ideal for the individual who has an active lifestyle, is "at low risk" for, or has a history of pressure ulcers.
- Available in dual and four valve/compartment to address seating challenges such as pelvic obliquities and sacral seating, to enhance postural control and to improve overall function.
- Light Weight, Approximately 3 pounds/1.36 kilograms.
- Available in Pediatric, Adult and Bariatric Sizes.
- 24 Month Manufacturers Limited Warranty.
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene.
- No Weight Limit