



*Solutions For Life*



POWERED BY  coolcore®

# GLACIAL SP



E2603/E2604

### Features:

- Contoured, high-density molded foam base
- Soft gel top layer for skin protection
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- 4-Way Stretch Poly or Urethane liner optional
- Sizes available - 10"-24" W x 10"-24"D
- Optional Mesh and Stealth Tek cover
- Optional ABS or ACM solid seat insert

### What is CoolCore®

- Chemical Free material • UPF/Flammability rated • 4-Way Stretch nylon interlock with Elastane
- Patented combination of chemical-free blends that deliver 3 distinct functions:  
**Moisture Wicking to Moisture Transportation to Moisture Evaporation.**
- **CoolCore® Technology is exclusive to Stealth Products.**



CHEMICAL-FREE THERMOREGULATION TECHNOLOGY

