



PREMIERE P



E2605/E2606

FEATURES:

- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base.
- Coccyx and sacral relief.
- FRG SunMate is an all natural flame retardant material.
- CoolCore® technology cover standard.
- Optional 4-Way Stretch Poly liner.
- Sizes available - 16"-22"W x 16"-22"D
- Standard Visco topper (exclusive FRG Sunmate) for added comfort
- Built-in positioning platform for off loading.
- Optional Urethane liner
- Optional mesh outer cover.
- Optional ABS or ACM Solid Seat insert.







PREMIERE P

STANDARD SIZES AVAILABLE

WIDTH	DEPTH			
	16" (40.64cm)	18" (45.72cm)	20" (50.8cm)	22" (55.88cm)
16"(40.64cm)	●	●	●	
18"(45.72cm)	●	●	●	
20"(50.8cm)	●	●	●	
22"(55.88cm)		Coming Soon	Coming Soon	Coming Soon

PRESSURE MAPPING



COOLCORE® TECHNOLOGY Exclusive to Stealth Products

- Chemical Free material • UPF/Flammability rated
- 4-Way Stretch nylon interlock with Elastane

Our patented combination of chemical-free blends deliver 3 distinct functions:
Moisture Wicking to Moisture Transportation to Moisture Evaporation.

