



Full Body System

Give your clients 50 hours of therapy a week – with TheraTogs!

The TheraTogs Full Body System was originally designed for children with complex neuromotor disorders, and is now available in a full range of sizes for a variety of applications.

TheraTogs are an effective modality for the management of neuromotor disorders caused by a wide range of diagnoses:

- Cerebral palsy
- Down Syndrome and other chromosomal abnormalities
- Traumatic Brain Injury – acute, post-acute, and chronic phases
- Hemiplegia due to Stroke – acute, post-acute, and chronic phases
- Multiple Sclerosis

The Full Body System allows the clinician to address the greatest variety of functional alignment and movement issues, extending from the torso and body core to the extremities proximal to the wrists and ankles.



If your client has these indications...	A Full Body System can provide...
Flexible kyphosis	Improved posture and respiratory capacity
Flexible lordosis / excess pelvic tilt	Improved posture
Difficulty maintaining midline trunk and sitting balance	Reduced deviations between thorax and pelvis; increased sensory awareness
Scissoring in stance or gait	Increased base of support, assist hip abductors
Balance problems due to hip muscle weakness	Assist to hip stabilizing musculature
Problems with proximal limb joint control, including elbow and knee	Reduced functional degrees of freedom; assist to muscles lacking adequate recruitment
Intoed or out-toed gait, or windblown gait	Reduced gait deviations
Recovery following orthopedic surgery at the spine, hip, or femur	Improved stability and comfort; strapping to support surgical goals



Full Body System



Model Sizes

Model No.	Description
#GFB100	Infant: Boys/girls 15-25 lbs
#GFB200	Preschool: Boys/girls 25-45 lbs
#GFB300	Pediatric: Boys/girls 45-65 lbs
#GFB400	Petite Adult M/F 65-110 lbs
#GFB600	Small Adult Female: 110-140 lbs
#GFB610	Medium Adult Female: 140-170 lbs
#GFB620	Large Adult Female: 170-200 lbs
#GFB700	Small Adult Male: 110-145 lbs
#GFB710	Medium Adult Male: 145-170 lbs
#GFB720	Large Adult Male: 170-190 lbs
#GFB730	X-Large Adult Male: 190-240 lbs

About TheraTogs

TheraTogs™ are FDA Class I medical devices intended to be issued by, and applied under the supervision of, a licensed healthcare practitioner engaged in neuromotor or sensorimotor retraining. TheraTogs systems are made of GoldTone™ - a proprietary and patented composite fabric with foam backing made of an aqueous-based elastomeric urethane. TogRite™ strapping consists of elastomeric strapping with an inert, silicone-based grip surface. TheraTogs components are not made with natural rubber latex.