



Defender Reha US

The Special Needs 2-in-1 CarSeat Solution
Helping Families Enjoy Safe Travelling

Thomashilfen & **harmony** are two family owned companies passionate about providing safe and comfortable products. We are committed to developing and delivering products that help families enjoy safely travelling together.

Let us **move** you!

“The Defender Reha allows children with special needs 360° of impact protection. The head, trunk and thighs supports, combined with our unique accessories, offer many positioning options. The seat allows for growth while remaining narrow enough for a second car seat as well as another child to fit in the back seat.”

Positioning-options INCLUSIVE



Ultra plush seat pad made with durable premium fabrics and extra foam for maximum comfort is removable for machine washing



Lateral trunk supports for an upright, stable sitting position

Thomashilfen

9-position, 1-hand height adjustable headrest

Thomashilfen



Head support holds the child's head safely in place

Deep, supportive headrest, backrest and seat base for 360° impact protection



Chest clip and soft harness covers included



Lateral side supports are removable to adjust seat width for plenty of room for the child to grow with the seat

The no-threaded 5-point harness provides additional support in Booster Mode for children weighing 30 to 110 lbs. Small crotch pad included.

Defender Reha US

The Positioning 2-in-1 Special Needs CarSeat

Cup/snack holder swivels and can easily be installed on either armrest



Smooth up-front harness adjuster

harmony

Accessories & Measurements



Footrest, short with footrest adapter and footrest pad



Footrest, long with footrest adapter



Seat wedge, below (+ 10°)



Table = playtray and additional impact shield



Abduction block



Crotch pad, large

Defender Reha **AND** all accessories exceed US FMVSS 213 and Canadian CMVSS 213 standards



Measurements & weights	Defender Reha US
Seat depth	13.5"
Seat width*	8.5 - 12"
Seat to head**	23.5 - 30.5"
Shoulder height (child)***	10 - 19.5"
Lower leg length	(7") 11 - 17"
Overall dimensions (w x h x d)	18.25 x 26.5 x 20"
Weight of seat	22 lb
HARNES MODE	
Weight of child	22 - 65 lb
Height of child	27 - 57"
BOOSTER MODE	
Weight of child	30 - 110 lb
Height of child	34 - 57"

* To reach max. seat width remove additional hip guide pads from additional seat pad
 ** To reach max. back height remove additional seat pad and foam pad under seat cushion
 *** To reach max. shoulder height remove additional seat pad and foam pad under seat cushion

Product & accessories	item code
Defender Reha US	886
Footrest, short (depth 4.7")* (tiltable and adjustable in angle)	806
Footrest pad (height 4") – (for 806)	834
Footrest, long (depth 7.1")* (tiltable and adjustable in angle)	805
Footrest adapter + Base plate*	821
Seat wedge, below (+ 10°), black	810
Table, black	831
Crotch pad, large, black**	832
Abduction block, black	833

* The footrest adapter must be ordered to use footrest
 ** We recommend using larger crotch pad only in booster mode.

