



# **RECARO** Monza Nova 2 Reha

Booster seat for children with special needs

”

*The safety of my child is most important to me. Generally, I am a careful driver, but you never know. In case of an accident I have to be 100% sure that the car seat protects my child against injuries.*

Susanne Bessel, mother



“

## **RECARO** Monza Nova 2 Reha is reliably tested many times

- ✓ FMVSS 213 (*NHTSA Standard for USA Car Seats*)
- ✓ UN/ECE R 44/04  
(*United Nations Economic Commission for Europe – regulation for child restraint systems*)
- ✓ DEKRA (*tests for „Special Needs Restraint“*)  
incl. testing of all rehab accessories
- ✓ Manufacturer’s guarantee



### **Pay attention to this seal**

Thus you can be sure that the **rehab seat** successfully passed the necessary overload-test up to 50 kg / 110.2 lbs according to ECE R 44/04.

**RECARO**  
**Monza Nova 2 Reha**  
 Rückhalteinrichtung  
 für spezielle  
 Anforderungen  
*restraint system  
 for special needs*  
 S  
 15 - 50 kg  
 33.1 - 110.2 lb

---

**CE**

- Head padding
- 5-point positioning harness
- Lateral trunk supports
- Crotch pad
- Lateral side supports



Rehab-equipment  
**INCLUSIVE** 

## ACCESSORIES



The practical padded swivel base, the Seatfix connection for ISOFIX and the footrest give additional stability.



The seat wedge inside and below provide a relaxed resting position during the journey. The abduction block offers extra positioning for children whose legs cross or "scissor".



The table is a support surface and an additional impact shield that also prevents the child from attempting to unbuckle the 5-point positioning harness.



**EXTRA SUPPORTIVE**  
 for children with special needs

# RECARO Monza Nova 2 Reha

Lightweight booster seat exceeds highest safety standard



Quick and easy seat transfer from one vehicle to another.



Improved strength and static-free side pads reduce crash-generated forces to the head, neck and shoulders by up to 30%.

**30%**



Lightest rehab seat in its class that you can carry with one finger ...  
Yet it carries a child up to 15 to 50 kg / 33 to 110 lbs.

Extra-firm armrests for even greater protection in the event of side impact.

**Thomashilfen** Finely designed to make life easier

Measurements & weights	RECARO Monza Nova 2 Reha
Seat depth	30 cm / 11.8" (43 cm / 16.9" with seat depth extension)
Seat width	front: 29 cm / 11.4" rear: 24 cm / 9.4" (20 cm / 7.9" with seat wedge, inside)
Back height	60 - 77 cm / 23.6 - 30.3"
Overall dimensions (w x h x d)	54 x 66 x 42 cm / 21.3 x 26 x 16.5"
Weight of seat approx.	6 kg / 13.2 lbs // 7 kg / 15.4 lbs (Seatfix model)
Carrying capacity	15 - 50 kg / 33.1 - 110.2 lbs
Children height	94 - 150 cm / 37 - 59"



Products & accessories	item code
RECARO Monza Nova 2 Reha	840
RECARO Monza Nova 2 Seatfix Reha (with Isofix)	841
Abduction block	804
Swivel base with footrest adapter** (incl. adapter plate)	820/1
<b>Recommended options for stability:</b>	
▪ Seatfix connection (for swivel base) <b>option 1</b>	818
▪ Stabilising bow, rear (for swivel base) <b>option 2</b>	819
Stand, front (for swivel base)	817
Footrest, long* (depth: 18 cm / 7.1") (tiltable and adjustable in angle)	805
Footrest, short* (depth: 12 cm / 4.7") (tiltable and adjustable in angle)	806
Seat depth extension* (+ 13 cm / + 5.1")	808
Seat wedge, inside	809
Seat wedge, below** (+ 10°)	810
Table	811
Soft harness strap cover, grey (pair)	815

\* Footrest and seat depth extension cannot be used in combination.

\*\* Swivel base and seat wedge (below) cannot be used in combination.