


ThevoSleepingStar A - The mobility bed for children

Thevo - The mobility bed for children

For a restorative sleep we need a balanced ratio between relaxation and movement. This is often disturbed with special needs children and leads to sleep disorders. With the mobility bed this deficit is being balanced.

ThevoSleepingStar A

- ✓ electronically driven actors provide the child with small movement impulses
- ✓ suitable for pain therapy, particularly for very immobile children
- ✓ suitable for children in vigil coma
- ✓ optimum pressure relief and movement support




Mattress

The mattress is designed especially to meet the needs of children. It is contaminant-free and as standard protected by an incontinence cover.

Control unit

Operates and controls the actors. A simple menu selection offers a variety of movement programmes. This provides for a stimulation individually adapted to the child.

Programmes:



ModelChoice

	Size	Item code
ThevoSleepingStar A incl. incontinence cover	70 x 140 cm / 27.6 x 55.1"	61311
ThevoSleeping Star A incl. incontinence cover	90 x 200 cm / 35.4 x 78.7"	65302/5
ThevoSleeping Star A incl. incontinence cover	100 x 200 cm / 39.4 x 78.7"	65322/5
Transporting and hygiene bag		64300
Bag for control unit		64302

Measurements & weights

	Size
Height of mattress	11 cm / 4.3"
Height of underframe	8 cm / 3.2"
Total height	19 cm / 7.5"
Carrying capacity	5 - 60 kg / 11 - 132.3 lb

ThevoSleepingStar cover



Incontinence cover

Features


- breathable
- cleaning by wiping disinfection
- sewn zipper
- washable at 65°C / 149°F
- suitable for tumble dryers
- autoclavable
- suitable for the re-use

Thevo
by Thomashilfen



For further information please contact your local dealer.

ThevoSleepingStar A



Thevo – The mobility bed

Thevo
by Thomashilfen

IGAP
Institute for Innovations in
Healthcare and
Applied Nursing Science



Item code: 95736, status: 09/2009, subject to technical alterations