

# EASyS Modular S – 2023 US Retail prices

|         |       |     |          |
|---------|-------|-----|----------|
| Date    | Buyer | PO# | Mark For |
| Address |       |     |          |
| Contact | Phone | Fax |          |

## Step 1 – Seat Sizes (Required options are indicated below by † symbol)

£ Modular S seat    £ Modular S seat



|                      | size 1 | size 2 | HCPCS  | Description |       |  |
|----------------------|--------|--------|--------|-------------|-------|--|
| †Seat                | 661011 | \$420  | 671011 | \$420       | E2292 | Solid contoured, integrated seat/back, 3 different head rest options for more head control; 140° back recline, growth adjustable; seat depth/ width adjustable from underneath, articulated foot plate height and angle adjustable; front/rear facing seat, easy to change direction; ventilated/washable upholstery ; knee angle adjustment with quick release lever; cozy comfortable extra seat cushion; calf panel; 5 different color options; tilt-in-space chassis available |
| †Back                | 661012 | \$518  | 671012 | \$518       | E2293 |  |
| †Recline             | 661013 | \$415  | 671013 | \$415       | E1014 |  |
| †Hip sup.            | 661014 | \$230  | 671014 | \$230       | E0956 |  |
| †Head rest bracket   | 661015 | \$170  | 671015 | \$170       | E1028 |  |
| †Calf panel          | 6635   | \$45   | 6735   | \$45        | E0995 |  |
| †Foot plate, flip up | 661016 | \$90   | 671016 | \$90        | K0040 |  |
| †Legrest, elevating  | 661017 | \$200  | 671017 | \$200       | E0990 |  |

| Seat Size Chart | Lower Leg   | Seat Depth  | Seat Width   | Back Height# | Max Height | Hip Angle  |
|-----------------|-------------|-------------|--------------|--------------|------------|------------|
| Size 1          | 7.3 – 12"   | 7.5 – 11.8" | 7.5 – 10.2"  | 11 – 16.5"   | 50"        | 90° - 140° |
| Size 2          | 7.3 – 14.6" | 10.6 – 15"  | 10.6 – 13.8" | 13.4 – 20.9" | 60"        | 90° - 140° |

  

| †Upholstery colors | Red/ Grey | Blue/ Grey | Green/ Grey | Magenta/ Grey | Grey/ Grey |
|--------------------|-----------|------------|-------------|---------------|------------|
| Size 1             | £ 6242    | £ 6243     | £ 6244      | £ 6241        | £ 6240     |
| Size 2             | £ 6442    | £ 6443     | £ 6444      | £ 6441        | £ 6440     |

†Head rests size 1&2 (Required option) HCPCS E0955

£ 62611 - Shell \$285

£ 62621 - Neck support \$285

£ 6267 - Angle adjust. laterals \$335    £ 6268 - fitting, ball joint \$90



## Step 2a – Outdoor Frame Option (Required options are indicated below by † symbol)

**Advantage-chassis** - Folding aluminum frame, tilt-in-space, back frame suspension, face-to-face option, 0° safety lock out, cargo basket, unfold lock



|                            | HCPCS | SIZE 1            | SIZE 2            |
|----------------------------|-------|-------------------|-------------------|
| Advantage-chassis          | E1234 | £ 62023    \$2295 | £ 64023    \$2295 |
| †Castor pin lock, pair     | K0073 | 67830    \$96     | 67830    \$96     |
| †Bus tie-down loop kit (*) | K0108 | 67839    \$205    | 67839    \$205    |
| Pneumatic wheels, full set | K0070 | £ 66461    \$279  | £ 66461    \$279  |
| Frame padding, black       | K0108 | £ 6622    \$125   | £ 6622    \$125   |
| Frame padding, red/ black  | K0108 | £ 66221    \$125  | £ 66221    \$125  |
| Tip assist/ curb tippers   | K0108 | £ 67502    \$99   | £ 67502    \$99   |
| Under seat storage basket  | K0108 | NC                | NC                |



| EASyS Modular S Technical Specs | SIZE 1              | SIZE 2       |
|---------------------------------|---------------------|--------------|
| Seat capacity                   | 77.2 lbs            | 88.2 lbs     |
| Frame capacity                  | 110.2 lbs           | 121.3 lbs    |
| Unit weight                     | 41 lbs              | 48 lbs       |
| Max seat tilt                   | 35°                 |              |
| Frame suspension                | Yes                 |              |
| Upholstery                      | Trevira-CS          |              |
| Foot angle                      | 78° - 100°          |              |
| Size of foot rest               | 11" x 7.1"          | 12.6" x 8.3" |
| Floor to handle                 | 28.2 – 46.1"        | 28.2 – 46.1" |
| Floor to seat                   | 19"                 | 19.5"        |
| Max width                       | 23.4"               | 27.2"        |
| Fold length                     | 33.5"               | 34"          |
| Fold height                     | 23.6"               | 27.2"        |
| Brakes                          | foot operated       |              |
| Seat direction                  | front & rear facing |              |
| Front wheels                    | 7.1" solid          |              |
| Rear wheels                     | 10" solid           |              |

\* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

## EASyS Modular S – 2023 US Retail prices

### Step 2b – Indoor Frame Option (Required options are indicated below by † symbol)

**Q-chassis** - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas springs, lockable rear wheels



|   | HCPCS | SIZE 1 & 2 |        |
|---|-------|------------|--------|
| Q-chassis indoor hi/low, pedal controlled | E1399 | £ 6666     | \$2198 |
| †EASyS Adapter for Q-chassis              | K0108 | £ 62951    | \$130  |
| Push bar for Q-chassis                    | K0108 | £ 6660     | \$329  |

#### Q-chassis Technical Specs

|                |              |
|----------------|--------------|
| Capacity       | 88.2 lbs     |
| Seat tilt      | -7 to +31°   |
| Floor to seat  | 11.4 - 27.3" |
| Seat direction | front facing |

### Step 3 – Positioning & General Options (Required options are indicated below by † symbol)

| Lower Body Positioning      | HCPCS | Unit | Small   | Medium   | Large   | Price |
|-----------------------------|-------|------|---------|----------|---------|-------|
| Foot straps with cam buckle | E0952 | Pair |         | £ 6648   |         | \$35  |
| Ankle Huggers, Bodypoint®   | K0108 | Pair |         | £ 213151 |         | \$205 |
| Footrest pad, height 1¾"    | K0108 | Each | £ 10622 |          | £ 10023 | \$99  |

| Mid Body Positioning                                       | HCPCS                            | Unit  | Small   | Medium | Large   | Price    |       |
|--|----------------------------------|-------|---------|--------|---------|----------|-------|
| Abduction block, padded                                    | E0957                            | Each  | £ 6620  |        | £ 6720  | \$160    |       |
| Hip supports   | 2 point, padded, SR buckle       | E0978 | Each    | £ 8420 |         | £ 8470   | \$95  |
|  | 4 point, padded, SR buckle       | E0978 | Each    | £ 8421 |         | £ 8471   | \$135 |
|  | Soft pelvic harness/ leg harness | E0978 | Each    | £ 6614 | £ 6714  | £ 213152 | \$199 |
| Crotch strap (fits 6629,8420,8470,8421,8471)               | E0978                            | Each  |         | £ 6341 |         | \$33     |       |
| Tray, height adjustable, clear with picture window         | E0950                            | Each  | £ 66281 |        | £ 67281 | \$190    |       |
| Mounting hardware tray - †required when ordering 6628/6728 | E1028                            | Each  | £ 66282 |        | £ 67282 | \$225    |       |

| Upper Body Positioning | HCPCS  | Unit  | Small | Medium   | Large    | Price  |       |
|------------------------|--|-------|-------|----------|----------|--------|-------|
| Chest harness supports | Butterfly style, neoprene                          | E0960 | Each  | £ 6385   |          | £ 6585 | \$199 |
|                        | Vest style   | E0960 | Each  | £ 6631   | £ 6632   | £ 6731 | \$199 |
|                        | H-style, posture support, small, Bodypoint®        | E0960 | Each  |          | £ 214806 |        | \$209 |
|                        | PivotFit shoulder harness, dynamic, Bodypoint®     | E0960 | Each  | £ 214270 | £ 214250 |        | \$209 |
| Trunk supports         | Rigid lateral, left & right, height adj.           | E0956 | Pair  | £ 6634   |          | £ 6734 | \$245 |
|                        | Flexible lateral l & r, w/chest strap, height adj. | E0956 | Pair  |          | £ 6633   |        | \$295 |

| General Accessories               | HCPCS  | Unit  | Small  | Medium | Large    | Price  |       |
|-----------------------------------|--|-------|--------|--------|----------|--------|-------|
| Grab bar, padded, black           | K0108  | Each  | £ 6621 |        | £ 6721   | \$140  |       |
| Footrest tube upholstery, grey    | K0108  | Pair  |        | £ 6230 |          | \$35   |       |
| Vent tray                         | E1029  | Each  | £ 6294 |        | £ 6494   | \$450  |       |
| Medical necessities pack, black   | K0108  | Each  |        | £ 6639 |          | \$80   |       |
| Oxygen tank holder, for D/E-tanks | E2208  | Each  |        | £ 6297 |          | \$210  |       |
| Straps                            | Secure folded seat/frame                     | K0108 | Each   |        | £ 667011 | \$30   |       |
|                                   | Wrist "runaway" tether                       | K0108 | Each   |        | £ 667111 | \$30   |       |
| Weather protection                | Canopy, sun/rain with window, black          | K0108 | Each   |        | £ 663761 | \$325  |       |
|                                   | Lower body windbreaker, black for 663761     | K0108 | Each   |        | £ 661852 | \$99   |       |
|                                   | Canopy, rain only, clear                     | K0108 | Each   |        | £ 6638   | \$245  |       |
|                                   | Rain shield, clear for canopy 663761/ 637811 | K0108 | Each   |        | £ 6739   | \$125  |       |
|                                   | Rain cape, grey                              | K0108 | Each   | £ 6340 |          | £ 6540 | \$163 |
|                                   | Sleeping bag, summer weight, black           | K0108 | Each   | £ 6619 |          | £ 6719 | \$299 |
|                                   | Sleeping bag, winter weight, black           | K0108 | Each   | £ 6618 |          | £ 6718 | \$299 |
| Umbrella, angle adjustable        | K0108  | Each  |        | £ 6594 |          | \$95   |       |

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.