



## Hercules Prime –2023 US Retail prices

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	


### Step 1 – Seat

Model	Price	HCPCS	Description
Booster seat	£ 461 \$3580	E1399	Booster seat with swivel (right vehicle door, incl. swivel lock) & tilt function (up to 10° adjustable), height adjustable backrest, height adj. head cushion, back angle adjustment for proper installation, seat depth adjustment, 5-point positioning harness with shoulder pads, crotch pad, mounting options for accessories, removable covers, machine washable, seat & accessories tested acc. FMVSS 213





Swivel function\* with lock out



Seat tilt up to 10°

This is a Universal Restraint System for disabled children, young persons and adults. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System. Only suitable for use in vehicles with 3-point (lap & shoulder) seat belts.

**This is a booster seat for children, young persons and adults weighing between 40 and 176.4 lbs. (18.1 and 80 kg) AND whose height is between 51.1 and 70.9 inches (130 cm to 180cm). Suitable only for forward-facing installation. THIS booster seat is NOT CERTIFIED for use in an aircraft.**




### Step 2 – † Required Options: **ALWAYS** use LATCH connectors (419) or stabilizing bow (421) for stability purposes!

	HCPCS	Unit	Item no.	Price
† Lower anchors connection	E1399	Each	£ 419	\$169
† Stabilizing bow (if car has no latch)	E1399	Each	£ 421	\$90



#### Accessories for additional stability purposes

Front stand	E1399	Each	£ 422	\$395
Cushion, leveling (for shaped seats)	E1399	Each	£ 513	\$150



#### Accessories for i ddYf VcXm

<YUXfYgřz \Y][ \hUX1 gtable, w. bracket	E1399	Each	£ 404	\$275
Shoulder guides, left & right, w. brackets	E1399	Pair	£ 427	\$245
Upper belt guide, additional*	E1399	Each	£ 425	\$61
Chest clip**	E1399	Each	£ 423	\$6



Technical Specs Hercules Prime	
Seat depth	40 - 46 cm / 15.7 - 18.1"
Depth behind crotch strap	24 - 38 cm / 9.5 - 15"
Seat width	40 cm / 15.7"
Seat to head	67 - 80 cm / 26.4 - 31.5"
Back width (seat)	45 cm / 17.7"
Shoulder height (belt guide)	42 - 64 cm / 16.5 - 25.2"
Shoulder width (with shoulder guide)	42 - 51 cm / 16.5 - 20.1"
Back angle (max.)	120°
Seat tilt-in-space (with tilting base 10°)	0° / 2.5° / 5° / 7.5° / 10°
Lower leg length (with foot rest)	17 - 43 cm / 6.7 - 16.9"
Height of lateral trunk supports (upper edge)	26 - 36 cm / 10.2 - 14.2"
Distance lateral trunk supports (max.)	35 cm / 13.8"
Max. overall dimensions (w x h x d)	52 x 88 x 43 cm / 20.5 x 34.6 x 16.9"
Min. folded dimensions (w x h x d)	52 x 21 x 73 cm / 20.5 x 8.3 x 28.7"
Weight of seat (basic model)	13.4 kg / 29.5 lbs
Body height user	130 - 180 cm / 51.1 - 70.9"
User weight	18.1 - 80 kg / 40 - 176.4 lbs

## Hercules Prime –2023 US Retail prices

### Accessories for trunk support

Lateral trunk supports, swing-away, w. brackets	E1399	Pair	£ 405	\$795
Lateral trunk supports, rigid, w. brackets	E1399	Pair	£ 429	\$245
Upper body guides, swing-away, w. brackets	E1399	Pair	£ 406	\$865
Upper body guides, rigid, w. brackets	E1399	Pair	£ 428	\$245



405 Lateral trunk supports, swing-away



429 Lateral trunk supports, rigid



406 Upper body guides, swing-away



428 Upper body guides, rigid

### Accessories for mid and lower body

Abduction block	E1399	Each	£ 412	\$223
Leg guides, left + right, w. brackets	E1399	Pair	£ 415	\$245
Footrest, foldable, angle adjustable	E1399	Each	£ 417	\$560



412 Abduction block



415 Leg guides, left + right, w. brackets



417 Footrest, foldable, angle adjustable

\* When using Hercules Prime on vehicle's **left back seat**, swivel function is **NOT** available!

#### \*\*CAUTION

Before installing chest clip you **MUST** determine if person needs chest clip installed **ABOVE** or **BELOW** strap adjusters. **ALWAYS** keep clip on chest height level. **Chest clip installations ABOVE chest height can cause serious injury!**

Hercules Prime **AND** accessories exceed US FMVSS 213 standards

Specs & pricing are subject to change without notice. HCPCS are suggestions only with no guarantee of reimbursement.  
ExoMotion, LLC dba Thomashilfen | Effective Date January 1st 2023