

Quote

 Order

Recaro Monza Nova 2 Latch Reha - 2019 US Retail prices

Date	Buyer	PO#	Mark For
Address			
Contact	Phone	Fax	



Seat w. Latch

Model	Price	HCPCS	Description
<input type="checkbox"/> 841-4	\$1299	E1399	Booster seat with Lower Anchors (additional anchor 818 needed to use swivel base combination), color carbon black booster seat, back with height-adjustable and inflatable headrest with integrated speakers, extra side impact protection, 5-point positioning belt, crotch pad, set of removable head, lateral trunk and side supports, seat and all accessories tested acc. FMVSS 213

Step 2 - General Options ("Either or" options are indicated below by † symbol)

General Accessories	HCPCS	Unit	Item no.	Price
Swivel base with footrest adapter**	E1399	Each	<input type="checkbox"/> 820-2	\$489
† Stabilizing bow**	E1399	Each	<input type="checkbox"/> 819	\$25
† Lower anchors connection**	E1399	Each	<input type="checkbox"/> 818	\$125
Front stand	E1399	Each	<input type="checkbox"/> 817	\$155
Footrest, size 2 (depth 7")* (tilt/ angle adjust.)	E1399	Each	<input type="checkbox"/> 805	\$270
Footrest, size 1 (depth 4.7")* (tilt/ angle adjust.)	E1399	Each	<input type="checkbox"/> 806	\$270
Seat depth extension (+ 5")*	E1399	Each	<input type="checkbox"/> 808	\$235
Seat wedge, inside	E1399	Each	<input type="checkbox"/> 809	\$135
Seat wedge, below (+10°)**	E1399	Each	<input type="checkbox"/> 810	\$170
Table	E1399	Each	<input type="checkbox"/> 811	\$280
Abduction block	E1399	Pair	<input type="checkbox"/> 804	\$105
Soft harness strap cover, black	E1399	Pair	<input type="checkbox"/> 815	\$37
Foot straps with cam buckle	E0952	Pair	<input type="checkbox"/> 6648	\$30

Technical Specs

Monza Reha	
Weight capacity	33.1 - 110.2 lbs
Seat weight	13.2 lbs
Upholstery	microfibre
Seat depth	11.8"
Seat width front	11.4"
Seat width rear	9.4"
Back height	23.6-30.3"
Dimensions	21.3x26x16.5"

This is a Universal Restraint System for disabled children. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System for this age group. Only suitable for use in vehicles with 3-point (lap & shoulder) seat belts. This is a Class II to III seat for children weighing between 33 and 110 lbs. (15 and 50 kg), and whose height is 59 inches (150 cm) or less. Suitable only for forward-facing installation.

* Footrest and seat depth extension cannot be used in combination

** Swivel base and seat wedge below or seat with lower anchors and seat wedge below cannot be used in combination

**** Swivel base requires either stabilizing bow 819 or lower anchors connection 818 for installation in a vehicle!**



Swivel base w. footrest adapter**



Footrest



Seat extension 5"



Seat wedge below +10°



Table



Abduction block and Seat wedge inside