

SWIFTY 2 – 2018 US Retail prices

| | | | |
|---------|----------------|--------------|----------|
| Date | Buyer | PO# | Mark For |
| Address | | | |
| Contact | Phone () | Fax () | |

Step 1 – Seat & Frame (Required options are indicated below by † symbol)



| | Swifty 2 | HCPCS | Description |
|---------------|----------|--------------|---|
| †Seat | 680821 | \$2825 E1236 | Folding frame with integrated/ contoured seat, fixed tilt, adjustable seat depth from underneath, footrest height adjustable, back recline, suitable for Early Intervention use |
| †Back | 680822 | \$125 E1014 | |
| †Bus tie-down | 6830 | \$175 K0108 | Bus tie-down loop kit * |
| H-belt | 6629 | \$5 E0960 | With purchase of Swifty |

| Seat Size Chart | Lower leg length | Seat depth | Seat width | Back height | Min. shoulder height | Hip angle |
|-----------------|------------------|------------|------------|-------------|----------------------|-------------|
| Swifty 2 | 8.3" – 14.6" | 14.8" | 15.6" | 27.6" | 11.8" | 90° to 100° |

Step 2 – Positioning & General Options

| Lower Body Positioning | HCPCS | Unit | Price |
|-----------------------------|-------|------|---------------------------------------|
| Foot straps with cam buckle | E0952 | Pair | <input type="checkbox"/> 6648 \$30 |
| Ankle Huggers** | K0108 | Pair | <input type="checkbox"/> 213151 \$205 |
| Footrest pad, height 1¾" | K0108 | Each | <input type="checkbox"/> 10023 \$95 |

Mid Body Positioning

| | | | | |
|--|----------------------------|-------|--------------------------------------|-------------------------------------|
| Abduction block, padded | E0957 | Each | <input type="checkbox"/> 6720 \$138 | |
| Hip supports | 2 point, padded, SR buckle | E0978 | Each | <input type="checkbox"/> 8470 \$71 |
| | Soft pelvic harness | E0978 | Each | <input type="checkbox"/> 6714 \$139 |
| Crotch strap (fits 6629,8470) | E0978 | Each | <input type="checkbox"/> 6341 \$25 | |
| Tray, height adjustable, clear with picture window | E0950 | Each | <input type="checkbox"/> 67281 \$190 | |
| Mounting hardware tray †required ordering 6728 | E1028 | Each | <input type="checkbox"/> 67282 \$205 | |

Upper Body Positioning

| | | | | |
|------------------------|--------------------------------------|-------|------|---------------------------------------|
| Chest harness supports | Butterfly style, neoprene | E0960 | Each | <input type="checkbox"/> 6585 \$145 |
| | Vest style | E0960 | Each | <input type="checkbox"/> 6632 \$145 |
| | H-style, w. integrated hip belt | E0960 | Each | <input type="checkbox"/> 66295 \$96 |
| | H-style, posture support, AEL | E0960 | Each | <input type="checkbox"/> 14813 \$185 |
| | PivotFit shoulder harness, dynamic** | E0960 | Each | <input type="checkbox"/> 214270 \$185 |
| | PivotFit shoulder harness** | E0960 | Each | <input type="checkbox"/> 214250 \$165 |
| Trunk supports | Rigid lateral | E0956 | Pair | <input type="checkbox"/> 6734 \$200 |
| | Flexible lateral with chest strap | E0956 | Pair | <input type="checkbox"/> 6633 \$285 |
| Head supports | Contoured pillow | E0955 | Each | <input type="checkbox"/> 6837 \$140 |

General Accessories

| | | | | |
|---|--------------------------------------|-------|-------------------------------------|---------------------------------------|
| Seat minimizer kit for Early Intervention use | K0108 | Set | <input type="checkbox"/> 6836 \$275 | |
| Pneumatic wheels, full set | K0070 | Set | <input type="checkbox"/> 6671 \$205 | |
| Cargo basket | K0108 | Each | <input type="checkbox"/> 6835 \$70 | |
| Frame padding, reversible, black/ grey | K0108 | Pair | <input type="checkbox"/> 6622 \$105 | |
| Grab bar, padded, black | K0108 | Each | <input type="checkbox"/> 6721 \$105 | |
| Tip assist | K0108 | Each | <input type="checkbox"/> 6829 \$69 | |
| Straps | Secure folded seat/frame | K0108 | Each | <input type="checkbox"/> 667011 \$25 |
| | Wrist "runaway" tether | K0108 | Each | <input type="checkbox"/> 667111 \$25 |
| Weather protection | Canopy, sun/rain with window, black | K0108 | Each | <input type="checkbox"/> 663761 \$295 |
| | Lower body windbreaker for 663761 | K0108 | Each | <input type="checkbox"/> 661852 \$95 |
| | Canopy, rain only, clear | K0108 | Each | <input type="checkbox"/> 6638 \$195 |
| | Rain shield, clear for canopy 663761 | K0108 | Each | <input type="checkbox"/> 6739 \$109 |
| | Rain cape, grey | K0108 | Each | <input type="checkbox"/> 6540 \$149 |
| | Sleeping bag, summer weight, black | K0108 | Each | <input type="checkbox"/> 6719 \$295 |
| | Sleeping bag, winter weight, black | K0108 | Each | <input type="checkbox"/> 6718 \$299 |
| Umbrella, angle adjustable, grey | K0108 | Each | <input type="checkbox"/> 6594 \$85 | |

Swifty 2 Technical Specs

| | |
|-----------------|---------------|
| Capacity | 110 lbs |
| Unit weight | 40 lbs |
| Seat tilt | 20° |
| Back recline | 90-100° |
| Upholstery | Trevira-CS |
| Floor to handle | 32.3"–48.8" |
| Floor to seat | |
| Max width | 26.4" |
| Fold length | 31.1" |
| Fold height | 19.3" |
| Brakes | foot operated |
| Seat direction | front facing |
| Front wheels | 7.1" solid |
| Rear wheels | 9.8" solid |



* Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.
 ** Bodypoint®