

## ThevoTherapy – 2018 US Retail prices

Date	Buyer	PO#	Mark For
Address			
Contact	Phone ( )	Fax ( )	

### Step 1 – Seat Sizes (Required options are indicated below by † symbol)



	<input type="checkbox"/> ThevoTherapy seat size 1	<input type="checkbox"/> ThevoTherapy seat size 2	HCPCS		
†Seat	840011 \$340	845011 \$380	E2294		Individual body adapted, integrated seat/back; 122° recline by gas spring; growth adjustable; articulated footrest height adjustable; ventilated/washable upholstery; front facing seat
†Back	840021 \$360	845021 \$420	E2293		
†Recline	840031 \$350	845031 \$350	E1014		
†Hip sup.	840032 \$200	845032 \$200	E0956		

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Seat to Head	Max Height	Hip Angle
Size 1	6 ¾ – 12"	7 ¼ – 12"	7 ¼ – 12 ¼"	18 ½ – 25 ¾"	54"	90° - 122°
Size 2	10 – 14 ¼"	10 ¾ – 15"	9 ¾ – 14 ¼"	25 ¼ – 30 ¾"	63"	90° - 122°

### Step 2 – Frame Options

**Hi-low chassis** - therapy frame, height/tilt easily adjustable by gas spring, 4 swiveling wheels, back wheels lockable

	HCPCS	Size 1	Size 2
Hi-low indoor chassis, pedal controlled gas spring	E1233 <input type="checkbox"/>	840041 \$1349 <input type="checkbox"/>	845041 \$1349 <input type="checkbox"/>
Push bar	K0108 <input type="checkbox"/>	8439 \$225 <input type="checkbox"/>	8439 \$225 <input type="checkbox"/>

### Step 3 – Positioning & General Options

Lower Body Positioning	HCPCS	Unit	Small >>>Medium >>>Large	Price
Foot straps with cam buckle	E0952	Pair	<input type="checkbox"/> 6648	\$30
Ankle Huggers, Bodypoint®	K0108	Pair	<input type="checkbox"/> 213151	\$205

#### Mid Body Positioning

Abduction block, padded	E0957	Each	<input type="checkbox"/> 6620 <input type="checkbox"/> 6720	\$138	
Hip supports	2 point, padded, SR buckle	E0978	Each	<input type="checkbox"/> 8420 <input type="checkbox"/> 8470	\$71
	4 point, padded, SR buckle	E0978	Each	<input type="checkbox"/> 8421 <input type="checkbox"/> 8471	\$170
Crotch pad (fits 8420,8470,8421,8471)	E0978	Each	<input type="checkbox"/> 6341	\$25	
Arm rest, height, depth and angle adjustable	E0973	Pair	<input type="checkbox"/> 8431 <input type="checkbox"/> 8478	\$375	
Tray, fold-away, black, height adjustable,	E0950	Each	<input type="checkbox"/> 8425 <input type="checkbox"/> 8475	\$315	
Hip support, high w/ upholstery	E0956	Pair	<input type="checkbox"/> 8423 <input type="checkbox"/> 8473	\$149	

#### Upper Body Positioning

Chest harness supports	Butterfly style, neoprene **	E0960	Each	<input type="checkbox"/> 6385 <input type="checkbox"/> 6585	\$145
	H-style, posture support, AEL **	E0960	Each	<input type="checkbox"/> 14813	\$185
	PivotFit shoulder harness, dynamic, Bodypoint® **	E0960	Each	<input type="checkbox"/> 214270	\$185
	PivotFit shoulder harness, Bodypoint® **	E0960	Each	<input type="checkbox"/> 214250	\$165
Trunk supports	Lateral, dynamic	E0956	Pair	<input type="checkbox"/> 8422 <input type="checkbox"/> 8472	\$255
	Rigid lateral	E0956	Pair	<input type="checkbox"/> 8426 <input type="checkbox"/> 8476	\$200
	Flexible lateral w/ chest strap	E0956	Pair	<input type="checkbox"/> 8427	\$285
Head support	Shell*	E0955	Each	<input type="checkbox"/> 84240	\$245
	Swing-away supp. bracket shell*	E1028	Each	<input type="checkbox"/> 84241	\$125
	Neck support*	E0955	Each	<input type="checkbox"/> 84290	\$245
	Swing-away supp. bracket neck*	E1028	Each	<input type="checkbox"/> 84291	\$125

#### General Accessories

Calf strap	K0108	Each	<input type="checkbox"/> 6635 <input type="checkbox"/> NC	\$35
Pelvis support	K0108	Each	<input type="checkbox"/> 8432 <input type="checkbox"/> 8479	\$89
Belt bow **	K0108	Each	<input type="checkbox"/> 6371 <input type="checkbox"/> 6571	\$119

\* Head rest with bracket is required option because seat back grows through head rest height adjustment

\*\* Belt bow is required for usage of chest harness supports

#### Technical Specs #

ThevoTherapy	Size 1	Size 2
Capacity	77 lbs	88 lbs
Unit weight	44 lbs	47 lbs
Max seat tilt	31°	31°
Seat direction	front facing	
Upholstery	Trevira-CS	
Floor to handle	31-41"	
Floor to seat	10-26"	
Max width	21"	
Max length	29"	
Fold width	21"	
Fold height	28"	
Fold length	29"	
Brakes	foot operated	
Front wheels	3.9"	
Rear wheels	3.9", lockable	

# More details in Technical Manual