

TiLite 2GX Order Form

| | | |
|------|---|---------|
| ZRW9 | Spinergy Design Lab LX Wheels (If same color is selected, anodized hubs and rims match color anodized package on page 1. Available with black spokes only.) <input type="checkbox"/> Gold Anodize <input type="checkbox"/> Green Anodize <input type="checkbox"/> Blue Anodize <input type="checkbox"/> Red Anodize <input type="checkbox"/> Purple Anodize <input type="checkbox"/> Orange Anodize..... | \$1,364 |
| | <input type="checkbox"/> 22"(501) - Std. Hub | |
| | <input type="checkbox"/> 24"(540) - Std. Hub | |
| | <input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.) | |
| | <input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.) | |

| | | |
|-------|--|---------|
| XRW16 | Spinergy Carbon Blade (Not available with Heavy Duty Frame.)..... | \$1,602 |
| | <input type="checkbox"/> 24"(540) | |
| | <input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth B.) | |

| | | |
|-------|--|-----|
| XRW17 | X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Sliver Anodized Handrims.)..... | NCO |
| | <input type="checkbox"/> 22"(501) | |
| | <input type="checkbox"/> 24"(540) | |

AXLES

| | | |
|-------|---|-------|
| XAXL1 | <input type="checkbox"/> Stainless Quick Release..... | STD |
| XAXL2 | <input type="checkbox"/> Titanium Quick Release..... | \$277 |
| XAXL3 | <input type="checkbox"/> Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.)..... | \$169 |

PNEUMATIC REAR TIRES

| | | |
|--------|---|-------|
| XRTR1 | <input type="checkbox"/> Treaded 22"(501), 24"(540), and 26"(590) (Not available with 25"(559) Rear Wheels.)..... | STD |
| XRTR2 | <input type="checkbox"/> Airless Inserts..... | \$98 |
| XRTR31 | <input type="checkbox"/> New MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559) and 26"(590)..... | STD |
| XRTR30 | <input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Primo Sentinel..... | \$136 |
| XRTR8 | <input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution..... | \$152 |
| XRTR32 | <input type="checkbox"/> 24"(540) Primo C-763 Knobby..... | \$104 |

E2213

SOLID REAR TIRES

| | | |
|--------|---|-------|
| XRTR22 | <input type="checkbox"/> 22"(501), 24"(540) and 25"(559) Primo Express (Not available with Spinergy Carbon Blade Rear Wheels.)..... | \$120 |
| XRTR23 | <input type="checkbox"/> 24"(540) and 25"(559) Primo Active..... | \$147 |

HANDRIMS

| | | |
|--------|--|-------|
| XRIM1 | <input type="checkbox"/> Aluminum – Silver Anodized..... | STD |
| XRIM21 | <input type="checkbox"/> Aluminum – Black Anodized..... | NCO |
| XRIM12 | <input type="checkbox"/> Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.)..... | \$509 |
| XRIM10 | <input type="checkbox"/> The Surge (All-In-One Oval with Gription™ Strip, Larger Profile)..... | \$433 |
| XRIM11 | <input type="checkbox"/> The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile)..... | \$466 |
| XRIM20 | <input type="checkbox"/> Q-Grip (Exclusive High-Friction, Anti-Slip Coating)..... | \$358 |

The Natural-Fit® (Original, Larger Profile) (Not available with 20"(451) Rear Wheels.)

| | | |
|--------|---|-------|
| TRIM14 | <input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)..... | \$428 |
| TRIM15 | <input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)..... | \$466 |
| TRIM16 | <input type="checkbox"/> Without Thumb Grip..... | \$390 |

The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)

| | | |
|--------|---|-------|
| TRIM17 | <input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)..... | \$440 |
| TRIM18 | <input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)..... | \$514 |
| TRIM19 | <input type="checkbox"/> Without Thumb Grip..... | \$423 |

Additional Handrim Options

| | | |
|--------|--|-------|
| XRIM5 | <input type="checkbox"/> Plastic Coated – Black..... | \$174 |
| XRIM22 | <input type="checkbox"/> Omit Handrims..... | NCO |

Tab Length

| | | |
|-------|---|------|
| XRIM7 | <input type="checkbox"/> Long Tabs..... | STD |
| XRIM8 | <input type="checkbox"/> Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.)..... | \$28 |

WHEEL LOCKS (Composite wheel locks not always compatible with 12"-13" Seat Depth. Uni-locks will be substituted at no charge.)

| | | |
|--------|---|------|
| XWLK7 | <input type="checkbox"/> Compact Composite Push to Lock..... | NCO |
| XWLK8 | <input type="checkbox"/> Compact Composite Pull to Lock..... | NCO |
| XWLK10 | <input type="checkbox"/> Extension Handles for Compact Composite Push/Pull to Lock..... | \$93 |

E0961

Additional Wheel Locks

| | | |
|--------|---|-------|
| XWLK11 | <input type="checkbox"/> MBL Composite Scissor Lock (Not available on 2GX Swing Away.)..... | \$158 |
| XWLK3 | Uni-Lock Set At: <input type="checkbox"/> Push OR <input type="checkbox"/> Pull..... | \$179 |
| XWLK4 | <input type="checkbox"/> Extension Handles for Uni-Lock..... | \$93 |
| XWLK15 | <input type="checkbox"/> Grade Aids Push to Lock (Not available with solid rear tires)..... | \$190 |

E0961

E0974

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

| | | |
|-------|---|-----|
| XBUP1 | <input type="checkbox"/> Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.)..... | STD |
|-------|---|-----|

TiLite 2GX Order Form

Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)

| | | | |
|-------|---|-------|-------|
| XBUP3 | <input type="checkbox"/> Padded Nylon..... | \$369 | E2611 |
| XBUP5 | <input type="checkbox"/> AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.)..... | \$369 | E2611 |
| XUPC3 | Ultrasuede Backrest Stripe <input type="checkbox"/> Black (Not available with Padded Nylon or AIR Breathable Back Upholstery.)..... | \$115 | |
| XBUP0 | <input type="checkbox"/> Omit Back Upholstery..... | NCO | |

SEAT UPHOLSTERY

| | | | |
|-------|---|-----|--|
| XSUP4 | <input type="checkbox"/> Tension Adjustable (Black only)..... | STD | |
| XSUP0 | <input type="checkbox"/> Omit Seat Upholstery..... | NCO | |

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

| | | | |
|--------|--|-------|-------|
| XCSH2 | <input type="checkbox"/> Comfort Curve Cushion (Stretch-Air cover)..... | \$122 | E2601 |
| ZCSH23 | <input type="checkbox"/> Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion..... | \$456 | E2607 |
| XCSH15 | <input type="checkbox"/> ROHO® Hybrid Elite SR™ Single Compartment Cushion (Fits seat widths 14" to 19")..... | \$552 | E2622 |
| XCSH14 | <input type="checkbox"/> ROHO® QUADTRO SELECT® Cushion (Fits seat widths 14" to 19") Profile: <input type="checkbox"/> Low <input type="checkbox"/> Mid <input type="checkbox"/> High..... | \$594 | E2624 |
| XCSH16 | <input type="checkbox"/> ROHO® Hybrid Select Cushion (Fits seat widths 14" to 18", and 20")..... | \$571 | E2624 |
| XCSH17 | <input type="checkbox"/> Vicair Active O2 Low (Fits even seat widths A 14" to 20". Non-stock upcharges may apply)..... | \$515 | E2624 |

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS

Aluminum

| | | | |
|-------|--|-------|--|
| XSDG2 | <input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std..... | \$217 | |
| XSDG3 | <input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std..... | \$233 | |
| XSDG6 | <input type="checkbox"/> Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.)..... | \$363 | |

2x2 Twill Weave Carbon Fiber

| | | | |
|-------|---|-------|--|
| XSDG7 | <input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std..... | \$336 | |
| XSDG8 | <input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std..... | \$347 | |

Custom Offset Profiled Sideguards CMD 007. Please specify desired offset to nearest ±1/2" (e.g., ±0.50") and sideguard type (rigid or fold-down) in the Notes section of the order form. Sideguard material must be selected. Some limitations may apply.

ARMRESTS (Desk Arms come with built-in Side Guards.)

| | | | |
|-------|--|-------|-------|
| XARM1 | Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.)..... | \$244 | E0973 |
| | <input type="checkbox"/> 9" - 12" OR <input type="checkbox"/> 12" - 15" | | |
| XARM4 | <input type="checkbox"/> Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.)..... | \$342 | E0973 |
| XARM2 | Removable Height Adjustable DeskArm with Rigid Side Guard (Requires 14" minimum Seat Depth B.)..... | \$363 | E0973 |
| | <input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15" | | |
| | <input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14") | | |
| XARM3 | Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth B.)..... | \$406 | E0973 |
| | <input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15" | | |
| | <input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14") | | |

Available Armrest Option

| | | | |
|-------|--|-------|--|
| XARM5 | <input type="checkbox"/> Covered in Black Ultrasuede (Not available with Desk Arms)..... | \$120 | |
|-------|--|-------|--|

SEAT BELTS

| | | | |
|--------|---|-------|-------|
| XBLT1 | <input type="checkbox"/> Velcro Adjustable, Black Nylon Webbing..... | \$61 | E0978 |
| XBLT2 | <input type="checkbox"/> Bodypoint® 2" Push-Button Buckle (Auto-Style.)..... | \$103 | E0978 |
| XBLT4 | <input type="checkbox"/> Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.)..... | \$196 | E0978 |
| XBLT5 | <input type="checkbox"/> Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.)..... | \$185 | E0978 |
| XBL10 | <input type="checkbox"/> Bodypoint® Optional Band Clamp Mounting Hardware..... | \$22 | K0108 |
| XBLT16 | <input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Push Button..... | \$201 | E0978 |

CALF STRAP

| | | | |
|-------|---|-------|-------|
| XCLF2 | Velcro Adjustable..... | \$44 | K0038 |
| | <input type="checkbox"/> Small (Recommended for 9" - 12" Footrest Widths) | | |
| | <input type="checkbox"/> Medium (Recommended for 11" - 15" Footrest Widths) | | |
| | <input type="checkbox"/> Large (Recommended for 14" - 17" Footrest Widths) | | |
| | <input type="checkbox"/> Xtra Large (Recommended for 16" - 19" Footrest Widths) | | |
| XCLF3 | Bodypoint® Padded Velcro Adjustable..... | \$109 | K0038 |
| | <input type="checkbox"/> Medium (Recommended for 12" - 15" Footrest Widths) | | |
| | <input type="checkbox"/> Large (Recommended for 16" - 20" Footrest Widths) | | |
| XCLF1 | <input type="checkbox"/> Buckle Adjustable, Black Nylon Webbing..... | \$34 | K0038 |

TiLite 2GX Order Form

RESIDUAL LIMB SUPPORT

| | | | |
|-------|--|-------|-------|
| XRSL1 | <input type="checkbox"/> Residual Limb Support (8" width x 10" length)..... | \$428 | E1020 |
| XRSL2 | <input type="checkbox"/> Residual Limb Support Pair (8" width x 10" length)..... | \$855 | E1020 |

BACKPACK

| | | | |
|-------|---|-------|--|
| XBCK3 | <input type="checkbox"/> Permobil Medical Necessities Backpack..... | \$103 | |
|-------|---|-------|--|

SEAT POUCH

| | | | |
|-------|---|-------|--|
| XPCH3 | <input type="checkbox"/> Permobil Medical Necessities Seat Pouch..... | \$66 | |
| XPCH4 | <input type="checkbox"/> Bodypoint® Mobility Bag..... | \$109 | |

CUSTOM CONFIGURATION OPTIONS

SPOKE GUARDS

| | | | |
|-------|--|-------|-------|
| XGRD1 | <input type="checkbox"/> Clear with TiMotion Logo (Not available with Mag Rear Wheels.)..... | \$147 | K0065 |
|-------|--|-------|-------|

ANTI-TIPS

| | | | |
|-------|---|-------|-------|
| XTIP1 | <input type="checkbox"/> Rear Aluminum - Black..... | \$174 | E0971 |
| XTIP2 | <input type="checkbox"/> User Friendly Flip-up..... | \$363 | E0971 |
| XTIP9 | <input type="checkbox"/> Tip Assist..... | \$115 | E0971 |

ACCESSORIES

| | | | |
|--------|--|-------|-------|
| XUEST1 | Upper Extremity Support Tray <input type="checkbox"/> 12" to 14" OR <input type="checkbox"/> 13" to 15" Seat Widths A (Requires Armrests.)..... | \$315 | E0950 |
| XLUG1 | <input type="checkbox"/> Luggage Carriers (Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.)..... | \$158 | |
| XMPCT2 | Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.) <input type="checkbox"/> Black..... | \$115 | |
| XCH1 | <input type="checkbox"/> Crutch Holder..... | \$244 | |
| XAFP1 | 4" Logo Patch <input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines..... | \$115 | |

GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Strain® is a registered trademark of Q'Strain Inc.

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR