


My chair is so fast
my friends have
to run to keep up.

— Wyatt



TiLite Pilot

Ready for take off



Specifications


Frame features	Carbon frame with 1" lower aluminum tubing	Seat back height	11" to 13" with standard back
Seat width	8" - 14"	Center of gravity adjustment	.5" - 3.5"
Seat depth	9" - 15"	Transit options	Transit Tie-Down, Transport
Front seat height	12" - 16.5"	Transport weight	11 lb*
Rear seat height	12" - 16"	Weight limit	95 lb transit, 120 lb standard
Front frame angle	Unlimited	Seat to footrest	3" - 14.5" without high mount parts

*Specifications may vary based upon configuration.
*Weight varies by configuration and does not include rear wheels.

Frame colors



Anodized colors



permobil

MS0057_RevA_Pilot

permobil




Ready to grow

The TiLite Pilot has built-in forward growth adjustment to grow with your child and maintain optimal rear wheel access. This gives kids the most efficient propulsion possible without needing to adjust or relearn their push mechanics. The TiLite Pilot lets kids cruise their world like never before.

Designed to move

The clean frame design features a seat pan and side panels that provide easy access to the adjustment bolts, making changes simple and intuitive. Unlike on other chairs, all of these adjustments can be done without compromising access to the rear wheels.



Cleared for landing

The ultra-lightweight design and carbon fiber frame panels make it 50% lighter and a breeze for kids to navigate all the adventures of their great big world. Quick release wheels and transit tie downs mean that bringing the pilot everywhere is also simple for mom and dad.