



TR

FORM & FUNCTION
AS ONE.



The sweeping curves of the titanium dual-tube frame are the perfect fusion of striking design and driven performance. Lighter, more efficient, and more elegant than any chair anywhere. As TiLite's flagship rigid chair, the TR is simply the best of the best.



TECHNICAL SPECIFICATIONS.

Frame Features	Titanium Dual-Tube Rigid TiFit Limited Adjustability
Frame Options	Black Edition SuperLite Edition Power Adaptable Frame Reinforced Frame
Transport Weight*	9.3 lb.
Weight Limit	265 lb.
Total Chair Width*	22"
Seat Width	12" - 20"
Front Seat Width (Seat Taper)	4 Options
Seat Depth	12" - 20"
Ergonomic Seat	7 Options
Front Seat Height	15" - 21"
Rear Seat Height	12.5" - 21"
Camber	0°, 2°, 4°, 6°, 8°, or 12°
Front Frame Angle	60°, 65°, 70°, 75°, 80°, 85°, or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	5"
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Folding Angle Adjustable Seat Back
Aesthetics	6 Titanium Frame Finishes 16 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

To learn more about TiLite options that can take your chair's performance to the next level, visit tilite.com/optimize.



TiLite | The Ultimate Ride™ | © 2015 TiLite, LLC.
800.545.2266 | tilite.com | MS0016_RevB_TR3
TiLite wheelchairs are shown with available options.