



# TWIST

CHILD  
PRODIGY.



Designed from the ground up for the growing child, the TWIST encourages exploration, independence, and, yes, playfulness. There are no parts or growth kits to buy. Instead, the TWIST itself adjusts as your child grows.



# TWIST

## TECHNICAL SPECIFICATIONS.

Frame Features	Aluminum Mono-Tube Rigid Growable Full Adjustability
Frame Options	Abduction Frame
Transport Weight*	12 lb.
Weight Limit	165 lb.
Total Chair Width*	19"
Seat Width	8" - 15" (2" built-in Seat Width Adj.)
Seat Depth	8" - 16" (3" built-in Seat Depth Adj.)
Front Seat Height	13" - 21.5" (up to 2.25" of built-in Front Seat Height Adj.)
Rear Seat Height	12" - 19.5" (up to 3.5" of built-in Rear Seat Height Adj.)
Camber	0°, 2°, 4°, 6°, 8°, or 12°
Front Frame Angle	70°, 80°, or 90°
Center of Gravity Adjustment	6"
Front Wheels	23 Options
Rear Wheels	35 Options
Handrims	15 Options
Convenience	Transit Tie-Down Option Folding Angle Adjustable Seat Back
Aesthetics	16 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
HCPCS Code	E1235
Warranty	Lifetime on Frame One Year on TiLite Components

\*Weights and widths based on 12" x 12" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weights. For reference, a pair of standard 22" rear wheels weighs 7 lb.

## GROWTH. BUILT IN.

Other chairs say they grow, but new parts or growth kits are required. With the TWIST, growth is built in. Sized from 8" x 8" up to 15" x 16", the TWIST offers 2" of growth in seat width and 3" of growth in seat depth.

## CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

To learn more about TiLite options that can take your chair's performance to the next level, visit [tilite.com/optimize](http://tilite.com/optimize).



TiLite | The Ultimate Ride™ | © 2015 TiLite, LLC.  
800.545.2266 | [tilite.com](http://tilite.com) | MS0032\_RevB\_TWIST  
TiLite wheelchairs are shown with available options.