

toppen helmet

- Machine Washable in warm water...line dry only
- Much cooler to wear than many other helmets
- Lightweight and durable, will not flake or tear
- Has a comfortable threaded buckle, not a bulky snap buckle
- Available in infant to young adult sizes
- Latex Free, blue nylon exterior fabric, red cotton interior fabric



Measuring

Measure circumference around the head, just above the ears



Order Information

18" (46 cm)	#106046
19" (48 cm)	#106048
20" (51 cm)	#106051
21" (54 cm)	#106054
22" (56 cm)	#106056
23" (58 cm)	#106058
24" (61 cm)	#106061
25" (65 cm)	#106065

Other sizes are available on request

Finally, a helmet that is durable, lightweight and easily washable!
