

Ascent Hospital Bed Mattress



A firm hospital bed mattress for easy transfers.

The Ascent is a firm hospital bed mattress made with state-of-the-art progressive laminated foams to provide a healthy skin environment. The edges of the mattress are extra firm to prevent users slipping off the mattress. The mattress is designed with sleep technologies that promote excellent comfort and healthy bedding solutions.

Mattress core is manufactured in an FDA approved facility that has decades of manufacturing experience to ensure the highest quality.

Key Features

- **Border edge protection to prevent slipping**
- **Zipper construction for easy cover replacement**
- **No flipping required**
- **Low profile: 6" tall**
- **Meets federal fire protection laws**
- **Made in USA**

Two Cover Options



Natural Bamboo Quilted Cover Elements

- Comfortable and fresh sleeping environment.
- Soft, luster, and silky feel.
- Natural properties: anti-bacterial, deodorizing, and breathable.
- Comprised of a Joma wool filling, which absorbs and causes moisture to evaporate, for maximum comfort.



Stretch Fabric Cover Elements

- A unique, waterproof cover that feels soft, just like fabric.
- Can be cleaned with household cleaners.
- Approved for FDA-grade cleaners.
- Will protect your mattress for many years.