

Supernal Recliner Bed System

All models available in Twin 80*, Full 80, Queen & Dual King.



Features	Supernal Recliner	Supernal Hi-Low	Supernal 5
Head & Foot	✓	✓	✓
Hi-Low	✗	✓	✓
Trendelenburg	✗	✗	✓
Reverse Trendelenburg	✗	✗	✓
Cardiac Chair	✗	✗	✓
Wall Hugging	✓	✓	✗
European-Style Head Tilt	✓	✓	✗
Twin 80 Capacity	400 lbs.	400 lbs.	300 lbs.
Full / Queen Capacity	500 lbs.	500 lbs.	300 lbs.
Standard Hand Control	Wireless	Wireless	Wireless
Special Hand Controls	✗	✗	✗
Massage	✓	✓	✓
Hoyer Compatible†	Low profile	w/ Casters‡	Low profile
Environmental Control Unit (ECU) Support	✗	✗	✗

*Headboard adds 1 ½” to length on all sizes

†Contact our office for details (800-475-8122)

‡5” Locking Casters enable use of Hoyer; net height increase of 6”

Supernal Technical Specs



Technical Specs	Supernal Recliner	Supernal Hi-Low	Supernal 5
Deck Height	11.5"	10.5"	12"
Range	N/A	10"	12"
Head Angle	65°	65°	60°
Foot Angle	35°	35°	35°
Trendelenburg Angle	N/A	N/A	10°
Reverse Trendelenburg Angle	N/A	N/A	10°
Cardiac Chair Angle	N/A	N/A	70°
Twin 80 Frame Weight [¶]	151 lbs.	263 lbs.	190 lbs.
Full 80 Frame Weight	199 lbs.	335 lbs.	235 lbs.
Queen Frame Weight	204 lbs.	340 lbs.	250 lbs.
Current	DC	DC	DC
Duty Cycle	10%	10%	10%

¶Uncrated, without mattress