

# Carrie® Potty Seat by Tumble Forms®... All Grown Up!

Accommodates adolescents to small adults



Little kids soon become big kids. That's why our popular Carrie® Potty Seat is available in three sizes; from adolescents to young adults. The above pictured small adult-size seat is specially designed to meet the toileting needs of users with heights ranging 58" to 68" (147-173 cm) and weighing up to 130 pounds (68 kg).

All seats have the same design features that make the Carrie Potty Seat the popular choice in positioning

equipment for toileting needs. The larger seat for adolescents and small adults provides secure positioning while encouraging toileting through increased flexion. The seat is shaped to position legs in slight flexion and abduction to encourage bladder and bowel functions. It can be positioned over 18" institutional toilets, many residential toilets or used as a free-standing commode—so it's ideal for clinical, school or home use.

**Other Features Include:**

- Two comfortable recline positions of 12° or 25°
- Lateral supports and tray help provide trunk stability
- Flexible elastomer splash guard/seat ring protects seat and is removable for easy cleaning
- Built-in anti-thrust seat helps maintain pelvic positioning
- Tumble Forms brand high-density foam has a washable, seamless coating that is urine and stain-resistant, nontoxic and latex free
- Four lockable casters on base facilitates safety with transfers

**Carrie® Potty Seat Includes:**

- Adjustable tray
- Height adjustable footrest
- Removable splash guard
- Adjustable Headrest
- Hip and chest harness
- 4 position removable footrest
- Removable bucket for use as a commode

**2771P** Elementary Carrie Potty Seat\*

**2771J** Junior Carrie Potty Seat\*

**2771SA** Small Adult Carrie Potty Seat\*

**Carrie® Potty Seat Specifications**

Item #	Size	Standing Height	Weight of Child
<b>2771P*</b>	Elementary	48" (125 cm)	60 lbs. (27 kg)
<b>2771J*</b>	Junior	58" (147 cm)	100 lbs. (45 kg)
<b>2771SA*</b>	Small Adult	68" (173 cm)	130 lbs. (68 kg)

\*Oversized and/or overweight shipping surcharge applies. Call for shipping quote.

*Note: may not fit over all toilets in reclined position.*



*Tray removes for ease in transferring*



*Splash guard removes for easy cleaning*



*Includes removable bucket*

**Available from your local Sammons Preston Rolyan  
Distributor or call Sammons Preston Rolyan**