


**Tumble
Forms®**

Vestibulator II[®] Flexi Disk



Flexi Disk Safety Notes

The Tumble Forms® 2772D Flexi Disk should only be used with the Tumble Forms® Vestibulator II® frame. The client weight should not exceed 150 lbs. Use floor mats. Never leave your client unattended. Please check your equipment regularly and replace any parts that show excessive wear. If there is any question of your clients safety, don't take a chance.

Flexion Swing

The Flexion Swing provides linear and orbital movement through space while the child sits on the platform and hugs the upright bolster portion of the swing.

