

# Carrie® Potty Seat by Tumble Forms 2®



Prices effective from January 1, 2014 through December 31, 2014.

Account No. \_\_\_\_\_  
Date: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Buyer: \_\_\_\_\_  
Attention: \_\_\_\_\_

Ship to: \_\_\_\_\_  
Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_  
Phone: (    ) \_\_\_\_\_

## Carrie® Potty Seat

- Tumble Forms 2 high-density foam has a washable, seamless coating that is urine and stain-resistant and nontoxic.
- Antimicrobial protection in the Tumble Forms 2 coating gives the seat an added level of protection against bacteria, mold and mildew, which cause stains, odor and product breakdown.
- Lateral supports and tray help provide trunk stability.
- Flexible elastomer splash guard/seat ring protects seat and is removable for easy cleaning.
- Built-in anti-thrust seat helps maintain pelvic positioning.
- Lockable caster base facilitates safety with transfers.

\_\_\_\_ **2771-P** Elementary      **\$2,588.00**  
\_\_\_\_ **2771-J** Junior          **\$3,014.60**  
\_\_\_\_ **2771-SA** Small Adult    **\$3,801.65**



Adjustable Tray  
& Footrest



Tray removes for ease  
in transferring



Splash guard removes  
for easy cleaning



Includes removable bucket

Lockable  
Casters



## Potty Seat Specifications

Size	Seat Width	Overall Width	Seat Depth	Overall Depth	Seat Height	Overall Height	Standing Height	Capacity
Elementary	10"	23"	10"	28"	20.5"	43.5"	48" (125 cm)	60 lbs. (27 kg)
Junior	12.5"	23"	13"	34.5"	21.5"	50"	58" (147 cm)	100 lbs. (45 kg)
Small Adult	14"	23"	15"	41"	32"	53.5"	68" (173 cm)	130 lbs. (68 kg)