

Tumble Forms 2[®] TriStander Order Form

Please mark clearly with an X the items you wish to purchase.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Features of the TriStander

- Adjustable height handle
- Light foot brake
- Adjustable, soft, cushioned pads
- Pads are seamless, impervious to body fluids and wipe clean
- Module cushions are interchangeable
- 3" double locking casters prevent roll and swivel
- Unique design allows prone, supine or vertical standing
- A foot-operated pneumatic tilt mechanism adjusts the stander from 15° to 90°
- Shoes adjust in width, height and angle
- Removable knee cushion

TriStander

___ 4525S	TriStander 45 (Yellow Frame)	\$3,961.70
___ 4530S	TriStander 58 (Blue Frame)	\$4,951.50
___ 4526S	TriStander 45 with Activity Tray	\$4,955.30
___ 4529S	TriStander 58 with Activity Tray	\$5,769.60

TriStander 45 Accessories

___ 4525F	Frame Only	21½"W x 38"L	\$2,809.05
___ 4525H	Headrest	9"W x 9¾"H	\$188.10
___ 4525K	Standard Trunk Support	8½" Inside W x 5½" Inside D	\$311.80
___ 4526K	Medium Trunk Support	10" Inside W x 6½" Inside D	\$334.55
___ 4525P	Hip Support	8½" Inside W x 4" Inside D	\$241.10
___ 4525N	Knee Supports (pr.)	3" Inside W x 2¼" Inside D	\$311.80
___ 4525E	Shoes with Bar	8½" Inside L x 2¾" H	\$575.65
___ 4525T	Activity Tray	23"W x 22"D	\$1,237.20
___ 9245-26	Head Support w/ Laterals	9"W x 9½"H, 3" side laterals	\$156.50
___ 9245-27	Adjustable Trunk Support	7½" Inside W x 6-7½" Inside D	\$486.05
___ 9245-28	Push Handle	7"L x 10"D	\$205.75

TriStander 58 Accessories

___ 4527F	Frame Only	23"W x 47"L	\$3,694.05
___ 4527H	Headrest	9"W x 10½"H	\$188.10
___ 4527K	Standard Trunk Support	11½" Inside W x 6" Inside D	\$352.30
___ 4526K	Medium Trunk Support	10" Inside W x 6½" Inside D	\$334.55
___ 4527P	Hip Support	12" Inside W x 6" Inside D	\$280.25
___ 4527N	Knee Supports (pr.)	4" Inside W x 2¼" Inside D	\$358.50
___ 4527E	Shoes with Bar	10½" Inside L x 3½" H	\$612.55
___ 4527T	Activity Tray	23"W x 22"D	\$1,465.75
___ 9245-26	Head Support w/ Laterals	9" W x 9½" H, 3" side laterals	\$156.50
___ 9245-27	Adjustable Trunk Support	7½" Inside W x 6-7½" Inside D	\$486.05
___ 9245-28	Push Handle	7"L x 10"D	\$205.75

Tumble Forms 2[®]
a Patterson Brand

With Anti-Microbial Protection
Odor Resistance, Increased Durability

