

HIPPOCAMPE

Beach and All-Terrain




Leisure activities everywhere and for everyone !

The Hippocampe Beach & All-Terrain accompanies you through all your outdoor activities. **Versatile and comfortable**, it adapts to your activities so you can enjoy moments of **relaxation and sharing with your relatives and friends**.

Thanks to its many options for comfort and leisure, the Hippocampe becomes a true companion to reach all terrains.

Enjoy hikes in the mountains or through forests, swims and walks in the snow, **discover new horizons !**



Characteristics

Sizes **S M L XL**

Colors **● ●**

Main leisure activities practiced

-Hike -Ski
-Swim -Surf -Travel
-Diving

▶ Customize your Hippocampe based on your specific needs

OPTIONS FOR YOUR COMFORT



Adjustable headrest

Adapts to the user's body type thanks to its height and depth adjustments.



Adjustable back

Adjustable in height and incline, it adapts perfectly to the user.



Adjustable armrests

Removable and rotating, they facilitate transfer to the chair.



Waist belt

Optimizes the position of the user in the wheelchair.



Mixed harness

In elastic neoprene, comfortable and water resistant, it facilitates the support and positioning of the user.



Anti-freezing cover

Comfortable and resistant it protects the user from cold and rain.



Protective nets

Adapt between the seat and the armrests or the footrest to protect the hands and feet of the user.



Parking brakes

Rotating and rustproof, they immobilize the wheelchair on an incline and facilitates transfers.

OPTIONS FOR YOUR LEISURE



Balloon wheels

Replace the double wheels and facilitate the wheelchair's mobility on sand for adult users.

Transport bag



The transport bag protects your wheelchair during your car and plane travels.

Front ski



Pair of back skis




Make your walks on snow easier, adapt to the terrain and prevent accumulation of snow on wheels.

Backrest bag



The backrest bag allows for the transport of your personal items during your leisure activities.



Seat length - feet

S : 72 cm - M : 82 cm
L : 92 cm - XL : 105 cm



Back seat height

37 cm to 47 cm



Seat width

40 cm



Seat depth

37 cm



Seat height - ground

34 cm




Interior width (without armrests)

50 cm



Total length

S : 170 cm - M : 180 cm
L : 190 cm - XL : 203 cm



Total height

103 cm



Total width

68 cm



Total length folded

S : 116 cm - M : 126 cm
L : 136 cm - XL : 149 cm



Weight

Weight without options : 17 kg
Max user weight : 130 kg



Size guide



Size S : to 1m40

Size M : from 1m40 to 1m65

Size L : from 1m65 to 1m85

Size XL : from 1m85

Zoom on the characteristics of the Hippocampe product range

Technical data

Seat : Comfortable hydrophobic foam (3cm thickness) covered by black fabric which dries quickly and avoids skin problems.

Frame : Stainless steel and anodized aluminum with a foam cover that is comfortable to the touch even in case of heat and cold exposure.

Push bar : Anodized aluminium tube with a handle covered in foam that improves comfort for the companion.

Certifications



Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA, ISO 13485, ISO 14971 and the NF santé.