

HIPPOCAMPE®

Pool 



Aquatic areas accessible to everyone

The Hippocampe Pool makes it easier for **people with reduced mobility to access aquatic and pool spaces.**

Usable in water, easy to use and hygienic **it accompanies users from changing rooms to pool areas** while avoiding multiple transfers.

The Hippocampe allows access to your pool and aquatic complex

Swimming is the most practiced sport by people with disabilities, it allows most individuals to use their motor skills in zero gravity situation.

Accessibility
Innovative Security
Rustproof Hygiene
Handling



Parking brakes included



Accessibility

The Hippocampe allows people with reduced mobility to move around freely between changing rooms, showers, and pool areas.

It can be submerged in water to simplify access to swimming.




Hygiene and security

The seat and footrest cover is made of comfortable hydrophobic foam that protects its user.


With easy maintenance, the Hippocampe can be rinsed and disinfected after each use.

FDA REGISTERED

Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA.




Seat length - feet
28.5" - 72 cm




Back seat height
14.5" - 37 cm



Seat width
15.5" - 40 cm



Seat depth
14.5" - 37 cm




Seat height - ground
13.5" - 34 cm



Interior width
19.5" - 50 cm



Total length
54" - 137 cm



Total height
38" - 96 cm



Total width
27" - 68 cm



Weight
Weight : 11.5 kg - 25,3 lb
Max user weight: 130 kg / 286 lb



Technical data

Seat : Comfortable hydrophobic foam (3cm thickness) covered by black fabric which dries quickly and avoids skin problems.

Frame : Stainless steel and anodized aluminum with a foam cover that is comfortable to the touch even in case of heat and cold exposure.

Push bar : Anodized aluminium tube with a handle covered in foam that improves comfort for the companion.

FDA REGISTERED

Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA.

