



CE 1000 BK



CE 1000 B



CE 1000 NBL






A



HCPCS
E1399

Adult Anterior Safety Walkers

Item #	Product Size	Color	Product Handlebar Height	Product Width Inside Handlebar	Product Dimensions (open at base)	Product Weight	Product Weight Capacity	UOM
CE 1000 B	Adult	Burgundy 	29" - 36"	17"	21.5" (W) x 25" (D)	16 lbs.	400 lbs.	1/bx
CE 1000 BK	Adult	Silver Vein 	29" - 36"	17"	21.5" (W) x 25" (D)	16 lbs.	400 lbs.	1/bx
CE 1000 NBL	Adult, Narrow Width	Blue 	28" - 35"	15"	20" (W) x 25" (D)	12 lbs.	400 lbs.	1/bx

- Height adjustable in 1" increments.
- Adjustable-brake spring tension.
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture.
- Skid-resistant, non-marring, easy-rolling neoprene wheels. Ideal for indoor and outdoor use.
- Folds easily and stands on its own in the folded position. (Figure A)

specifications

Carton Shipping Weight (CE 1000 B)	21 lbs.
Carton Shipping Weight (CE 1000 BK)	21 lbs.
Carton Shipping Weight (CE 1000 NBL)	17 lbs.
Warranty	Limited Lifetime