

# Exclusive VA Options and Accessories




# Custom Manual Options & ACCESSORIES

The **Quickie Service Appreciation Program** was designed to give back to those who have given so much to our country. Quickie has extended these exclusive features and options to our soldiers to express our gratitude and respect for the lives they have dedicated to the preservation of freedom and justice. Our family and friends sleep well at night and enjoy the liberty of the greatest country on earth because of your selfless acts.



**Show your military pride!** Six military service patches are available on the Jay J3 backrest and tension adjustable upholstery.



**Customize your look!**  
Choose from our exclusive to Quickie, Kolorfusion® patterns.

**Other camo options available through Built-4-Me**



### Quickie Xtender

- Power assist device to aid in propulsion
- No programming needed
- Easy to transport and store



### Sunparts Online

Quickly and accurately search and order parts for all Quickie® Zippie® and Breezy® wheelchairs, as well as JAY® seating products

Go to [www.sunrisemedical.com](http://www.sunrisemedical.com) and select the Parts tab to begin your search. Reverse lookup is available to find parts by entering the serial number.



COMPLEMENT YOUR QUICKIE OR ZIPPIE



# Casters and Forks

## CASTERS

### Aluminum Soft Roll



4" X 1.5"    5" X 1.5"    6" X 1.5"

### Aluminum Polyurethane / Soft Roll



4" X 1"    5" X 1"    5" X 1.5" Soft Roll

### Pneumatic




6"    8"    8" X 2" (Pneumatic with Insert)

### Semi-Pneumatic



4" X 1.25"    5" X 1.5"    6" X 1.5"    8" X 1.5"

### Low Profile Polyurethane




4"    5"    7"

### Polyurethane



6"    8"

### Micro and Micro Lighted



3" (Non-Lighted)    4" Lighted    5" Lighted

## FORKS



Multi Position Fork (Short pictured)



Standard Fork




Slotted Fork



Fork with Integrated Bearing 1 1/2" Stem Bolt (Tall Wide pictured) (Silver or Black)




3.5" Standard Fork



Standard Stem Bolt (on multi position fork)



3/4" Stem Bolt (on multi position fork)



1 1/2" Stem Bolt (on multi position fork)



Single Side Fork



Frog Leg



Caster Pin Lock

# Accessories

## POSITIONING BELTS



1.5" Auto Buckle Universal




2" Aircraft Buckle Positioning Belt Standard



2" Padded Aircraft Buckle Positioning Belt Padded

## SPOKE GUARDS



Quickie "Q" Clear



Quickie Black



Zippie Color

## BACK UPHOLSTERY OPTIONS



Standard Back



Tension Adjustable Back



Ballistic Back Upholstery



3DX™ Upholstery

NOTE: Not all accessories are available on all chairs.

## MISCELLANEOUS



Quickie Backpack



Kids Backpack




Quickie Seat Pouch



Padded Swing Away Armrest Detachable Gel Armpad Cover



Quickie Caddy




"L" Shaped Wrench Tool Kit

# Tires, Wheels and Handrims

## TIRES



## WHEELS



## HANDRIMS



## AXLES & VALVES




\*Available in Black, Blue, Red, Yellow, and White Spokes through Built-4-Me \*\*Available with Standard or Super Kevlar; Schwalbe, Spinerly, Natural Fit and Schrader are all registered trade names.


# Armrest and Hangers

## ARMRESTS Armrests shown in desk length



## HANGERS





# Footplates

# Wheel Locks and Side Guards

## FOOTPLATES - FOLDING



Composite



New Composite



New Angle Adjustable Composite



Adult Angle Adjustable



Kids Angle Adjustable



Locking Angle Adjustable



Aluminum Angle Adjustable



Split Platform Angle Adjustable Flip Away



Platform Flip Up



Tubular Flip Up (Available w/ Cover)



Auto Folding

## FOOTPLATES - RIGID



Tubular



Tubular with Flat Cover



Removable Snap On Cover



Aluminum Platform



Ultra Light Platform w/ Positioning Plates

## FOOTREST OPTIONS



Neoprene Impact Guards



Velcro Gel Cover for Hanger/Frame



Wear Pads



Positioning plates

## WHEEL LOCKS



High Mount Push Open



High Mount Push Closed



High Mount Pull Open



High Mount Pull Open



New High Mount Push Open Curved Handle



New High Mount Push Closed Curved Handle



New High Pull Open Curved Handle



New High Pull Closed Curved Handle



Ergo Scissor Open



Ergo Scissor Closed



Compact Open



Compact Closed




Rear Foot Operated Lock



Grade Aid



SR45 Extension Handles



Wheel Lock Extension Handles

## SIDE GUARDS



Plastic



New Plastic



Aluminum



Aluminum with Fender



Small Carbon Fiber