

zippie® Voyage®

GO TOGETHER,
GROW TOGETHER

early intervention device combines Baby Jogger™ stroller expertise, Zippie's incredibly versatile seating for growing babies, and one of the largest selections of functional accessories in the industry!



Check out the ZippieVoyage video

FEATURES AND BENEFITS

Baby Jogger City Select® Base

Intuitive, easy-to-use Baby Jogger City Select stroller features patented quick-fold technology, puncture-resistant rear wheels, and swivel front wheels with suspension and lock-out for long-distance strolls.



Carry-all Basket

Store up to 15 lbs. of your essentials below the seat.



Telescoping Stroller Handle

Allows for one-handed steering and adjusts to different caregiver heights.



NEW Slide-N-Lock™ Interface

Enables you to quickly and comfortably attach, remove, and reverse the seat in one simple motion.



Reversible Seat

Allows you to position the seat toward you to interact with your child, or forward so your child can explore their surroundings.



Optional Tilt & Recline

Assist with feeding, digestion, respiratory function, and visual orientation, and encourage achievement of developmental milestones such as bringing hands to the mouth, reaching for objects, or engaging in activities.



Options for Caregiver Convenience



Voyage Moderate Seating



Voyage Advanced Seating

Voyage Moderate Seating (VMS)

Designed for children with mild to moderate positioning needs. Offers soft, adjustable components and polka dot upholstery with colored accents.

Voyage Advanced Seating (VAS)

Designed for children with moderate to advanced positioning needs. Offers an extensive selection of seat options and adjustable components, and black stretch upholstery with colored piping.

TECHNICAL SPECIFICATIONS

HCPCS Code:

E1232
E1234 (without seating)

Frame Material:

Aluminum

User Weight Capacity:

75 lbs.

Transit Approved:

*Yes ISO 7176-19

FRAME AND FEATURE MEASUREMENTS

(Dependent on Configuration)

Overall Width:

24"

Seat Widths:

Advanced: 6" to 13"
Moderate: 9.5" to 12"

Seat Depths:

Advanced: 5" to 13"
Moderate: 8" to 13"

Back Height:

Advanced: 7" to 17"
Moderate: 16" to 24"
(includes head support)

Folded Height:

12"

Chair Weight:

Base only: 23 lbs.

Tilt Range:

35° (+5 to 30°)

Caster Wheels:

7"

Rear Wheels:

12" Pneumatic with Airless Insert

*Transit Approved:

Conforms with ISO 7176-19
for use in a motor vehicle